



promed®

Latex Powdered

Promed® Examination Gloves

- » Natural Rubber Latex
- » Non-Sterile
- » Single-Use
- » Powdered
- » Smooth Fingers

#### Comfort

*Natural Latex Rubber:*

- Natural elasticity provides a soft, comfortable fit and reduces hand fatigue; provides optimal sensitivity

*Lightly powdered:*

- USP grade cornstarch reduces sensitivity risk

*Ambidextrous:*

- Greater comfort, fit and feel

*Beaded Cuff:*

- Ease of donning and durability during use

*Cuff Length:*

- Available in standard 240mm

#### Protection

*Suitable Use:*

- Medical examination - non-invasive diagnostic and therapeutic procedures

*Barrier Protection:*

- Effective against potentially infective materials and other contaminants

*Strength and Durability:*

- High level protection; resistant to punctures with re-seal qualities

*Quality Standards:*

- CE Mark - Medical 93/42/EEC, ASTM3578-00, Exceeds AS4011 for pinholes, AQL1.5, ARTG No 109653

#### Latex Powdered

Promed® examination gloves

##### Product

Latex Powdered examination gloves; 100 gloves

##### Size

Small

##### Catalogue No.

PMD13000

##### Packaging

10/Carton

Latex Powdered examination gloves; 100 gloves

Medium

PMD13001

10/Carton

Latex Powdered examination gloves; 100 gloves

Large

PMD13002

10/Carton

Latex Powdered examination gloves; 100 gloves

X-Large

PMD13004

10/Carton

#### Glove Holder

One, two and three tier

##### Product

Glove holder, wall mounted

##### Size

One tier

##### Catalogue No.

APD3210

##### Packaging

Each

Glove holder, wall mounted

Two tier

APD3450

Each

Glove holder, wall mounted

Three tier

APD3640

Each



#### Glove Box Colours

##### Size

Small

##### Colour

Light Blue

Medium

Pink

Large

Blue

X-Large

Purple

ThermoFisher  
SCIENTIFIC

For customer service, call 1300-735-297  
To fax an order, use 1300-731-771  
Visit us online: [www.thermofisher.com.au](http://www.thermofisher.com.au)

©2010 Thermo Fisher Scientific Inc. All rights reserved. A.B.N. 52 058 390 917