



promed.

Nitrile Soft

Promed® Examination Gloves

- » Latex Free
- » Non-Sterile
- » Single-Use
- » Powder Free
- » 200 gloves/box for efficient storage, less waste and fewer dispenser changes
- » Blue coloured glove for easy identification

Comfort

Latex Free and Powder Free:

Improved Conformability:

Micro Textured Surface:

Ambidextrous:

Beaded Cuff:

Cuff Length:

- Reduces risk of allergic reactions and prevents type 1 allergic reactions
- Similar flexibility to natural rubber, conforms to the wearer's hand during use and allows optimal sensitivity.
- Provides non-slip grip in wet and dry conditions
- Greater comfort, fit and feel
- Ease of donning and durability during use
- Available in standard 240mm

Protection

Suitable Use:

Barrier Protection:

Strength and Durability:

Quality Standards:

- Medical examination - non-invasive diagnostic and therapeutic procedures
- Effective against potentially infective materials and other contaminants
- Resistant to puncture and tears; resistant to development of holes
- CE Mark - Medical 93/42/EEC, ASTM D 6319, Exceeds AS4011 for pinholes, AQL1.5, ARTG No 100975

Nitrile Soft

Promed® examination gloves

Product

Nitrile Soft examination gloves; 200 gloves
Nitrile Soft examination gloves; 200 gloves
Nitrile Soft examination gloves; 200 gloves
Nitrile Soft examination gloves; 200 gloves
Nitrile Soft examination gloves; 200 gloves

Size

Catalogue No.

Packaging

X-small	PMD173060	10/Cartron
Small	PMD173070	10/Cartron
Medium	PMD173080	10/Cartron
Large	PMD173090	10/Cartron
X-large	PMD173095	10/Cartron

Glove Holder

One, two and three tier

Product

Glove holder, wall mounted
Glove holder, wall mounted
Glove holder, wall mounted

Size

Catalogue No.

Packaging

One tier	APD3210	Each
Two tier	APD3450	Each
Three tier	APD3640	Each



Glove Box Colours

Size

Colour

X-Small	
Small	
Medium	
Large	
X-Large	

ThermoFisher
SCIENTIFIC

For customer service, call 1300-735-297
To fax an order, use 1300-731-771
Visit us online: www.thermofisher.com.au

©2010 Thermo Fisher Scientific Inc. All rights reserved. A.B.N. 52 058 390 917