



promed.

Vinyl Blue

Promed® Examination Gloves

- » Latex Free
- » Non-Sterile
- » Single-Use
- » Powder Free
- » Blue coloured glove preferred colour in food preparation

Comfort

Latex Free and Powder Free:

Beaded Cuff:

Cuff Length:

Ambidextrous:

Cuff Length:

- Reduces risk of allergic reactions and prevents type 1 allergic reactions
- Ease of donning and durability during use
- Available in standard 240mm
- Greater comfort, fit and feel
- Available in standard 240mm

Protection

Suitable Use:

High Chemical Resistance:

Anti-static Properties

Quality Standards:

- Medical examination where latex can not be tolerated
- Ideal for cleaning, catering, food handling and preparation
- Ability to withstand acids, fats and alcohol
- Reduces static
- CE Mark - Medical 93/42/EEC, ASTM D5250(01), ARTG No 150846

Vinyl Blue

Promed® examination gloves

Product

Vinyl Blue examination gloves; 100 gloves

Vinyl Blue examination gloves; 100 gloves

Vinyl Blue examination gloves; 100 gloves

Vinyl Blue examination gloves; 100 gloves

Size	Catalogue No.	Packaging
Small	PMD1200	10/Cartron
Medium	PMD1201	10/Cartron
Large	PMD1202	10/Cartron
X-large	PMD1203	10/Cartron

Glove Holder

One, two and three tier

Product

Glove holder, wall mounted

Glove holder, wall mounted

Glove holder, wall mounted

Size	Catalogue No.	Packaging
One tier	APD3210	Each
Two tier	APD3450	Each
Three tier	APD3640	Each

Glove Box Colours

Size	Colour
Small	
Medium	
Large	
X-Large	



ThermoFisher
SCIENTIFIC

For customer service, call 1300-735-297
To fax an order, use 1300-731-771
Visit us online: www.thermofisher.com.au

©2010 Thermo Fisher Scientific Inc. All rights reserved. A.B.N. 52 058 390 917