

## Housekeeping Products

### Natureclean

2



#### Neutral Hard Surface Cleaner

Natureclean is a revolutionary pH neutral product designed to be kind on the environment whilst effective in use. Natureclean contains ingredients which safely and effectively penetrate and remove soils without need for solvents or high alkalinity. Natureclean does not contain phosphates, silicates, formaldehyde or any form of alkali. Natureclean may be used for cleaning walls, floors and work areas. (This product is not intended for use on therapeutic devices).

**PACK SIZE: 2X5L P/C 17826**

**PACK SIZE: 15L P/C 16522**

### K-San

23



#### Disinfectant

K-San is a powerful hospital-grade disinfectant that does not require rinsing. K-San disinfects floors, walls, benches and processing equipment.

(This product is not intended for use on therapeutic devices).

**PACK SIZE: 2X5L P/C 17464**

### Easy

12



#### Creme Cleanser

Easy has been specially formulated to clean a wide range of surfaces without the harsh scratching some scouring powders can cause.

Easy is ideal for basins, tiles and stainless steel sinks and surrounds.

(This product is not intended for use on therapeutic devices).

**PACK SIZE: 2X5L P/C 16542**

**PACK SIZE: 12X1KG P/C 16540**

### Steelcoat

16



#### Stainless Steel Polish

Steelcoat has been specially formulated to protect and enhance stainless steel surfaces. Steelcoat is ideal for use in a range of applications where a high sheen, protective finish is required.

(This product is not intended for use on therapeutic devices).

**PACK SIZE: 2X5L P/C 16500**

## Miscellaneous Products

### Dermasoft



#### Antibac Hand Cleanser

Dermasoft Antibac Hand Cleanser is a hand cleaner containing antibacterial properties. Antibac is made under strict hygiene conditions for your assurance of quality and consistency. AUSTR 75220

**PACK SIZE: 8X1L P/C 6782801**

### Dermasoft



#### Sanitising Hand Gel

Dermasoft Sanitising Hand Gel is a waterless hand cleaning gel formulated to kill 99.99% (log 4 reduction) of organisms in 30 seconds including E.coli, MRSA, Listeria monocytogenes and Salmonella typhimurium. Sanitising Hand Gel contains emollients that protect the skin from drying out. AUSTR 61125

**PACK SIZE: 4X1L P/C 16458**

**PACK SIZE: 8X500mL P/C 17536**

**PACK SIZE: 12X100mL P/C 7354991**

### Dermasoft



#### Moisturiser and Barrier Cream

Dermasoft Moisturiser and Barrier Cream is a delicately fragranced hand moisturiser containing a light barrier cream. It quickly penetrates the skin leaving hands feeling soft and supple. Dermasoft Moisturiser and Barrier Cream leaves little residue on the hands and will not mark stainless steel.

**PACK SIZE: 4X1L P/C 17544**

## Accessories

### CSSD Test Kit

P/C 7370990



### Cap Spanner

63cm spanner aids in opening the wide neck 5L bottles. P/C 7380990



A Campbell Brothers Limited Business

cleaned to perfection


**ThermoFisher**  
SCIENTIFIC

For customer service, call 1300-735-297  
To fax an order, use 1300-731-771  
Visit us online: [www.thermofisher.com.au](http://www.thermofisher.com.au)

©2010 Thermo Fisher Scientific Inc. All rights reserved. A.B.N. 52 058 390 917

# Medical Cleaning System

Offering the complete cleaning solution  
for your Sterilising Facility



cleaned to perfection



# The Cleantec Medical Cleaning System

The Medical Cleaning Range offers exceptional quality and performance in sterilising facilities. The range provides a series of products that are efficient and easy to use, utilising concentrated products to improve economy.

Our whole range complies to AS/NZS 4187:2003 and AS/NZS 4185:2001.

## The Cleantec Medical Cleaning System features and benefits to you:

### A safer colour-coded system

All products are colour-coded and iconed to simplify identification for staff. Dispensing equipment has been designed to reduce handling.

Labels, Material Safety Data Sheets (MSDS) and wall charts comply with the latest Australian Standards.

### Economical and effective

The Medical Cleaning Range offers excellent value in terms of price and performance. Concentrated products and automatic dispensing gives the correct dosage to ensure consistent results.

### Product information and training

The correct use of each product is detailed on the colour fact sheets, customised wall charts and bottle labels. Our technical representatives provide training on how to optimise cleaning and hygiene and how to safely handle cleaning products. Material Safety Data Sheets are supplied for each product. The MSDS's are also available from our website at [www.cleantec.com.au](http://www.cleantec.com.au) for your convenience.



## Manual Cleaning

### Sonic

1



#### Ultra Sonic Detergent

Sonic is a mildly alkaline medical instrument and equipment detergent formulated to improve the cleaning rate and efficiency of ultrasonic cleaners. Sonic can also be used in pre-wash cleaning and manual cleaning of instruments and medical equipment and complies with AS/NZS 4187:2003 and AS/NZS 4185:2001.

Sonic is available either in a wide-neck or encapsulated, colour-coded and numbered bottle.

**PACK SIZE: 2X5L(WIDE-NECK) P/C 7322205**

**PACK SIZE: 2X5L P/C 17022**

**PACK SIZE: 15L P/C 7181015**

### Release Plus

1



#### Proteolytic Enzymatic Detergent

Release Plus is a concentrated multi-enzyme detergent formulated for manual cleaning and soaking of medical instruments and equipment.

Release Plus effectively removes blood, proteinaceous soils and lipids from glass, plastic and soft metal and complies with AS/NZS 4187:2003 and AS/NZS 4185:2001.

Release Plus is available in an iconed and numbered bottle.

**PACK SIZE: 2X5L P/C 6944205**

## Automatic Cleaning

### Achieve

3



#### Washer/Decontaminator Detergent

Achieve is a revolutionary instrument processing detergent developed to efficiently clean instruments, containers and articles, in a range of mechanical washers. In order to comply with the Australian Standards this product is designed to buffer at lower pH's than other washer detergent products. Achieve's advanced formula complies with AS/NZS 4187:2003 and AS/NZS 4185:2001.

Whilst being non-dangerous, Achieve is available in an encapsulated colour-coded, iconed and numbered bottle.

**PACK SIZE: 2X5L P/C 7344205**

### NewGen

3



#### Low Alkaline Decontamination Detergent

NewGen is a revolutionary instrument processing detergent developed to efficiently clean instruments, containers and articles, in a range of mechanical washers. In order to comply with the Australian Standards this product is designed to buffer at lower pH's than other washer detergent products. NewGen is suitable for washers and decontaminators with proprietary dispensing systems and complies with AS/NZS 4187:2003 and AS/NZS 4185:2001.

NewGen is the same formulation as Achieve, available in a wide-neck, colour-coded, iconed and numbered bottle.

**PACK SIZE: 2X5L P/C 7343205**

**PACK SIZE: 15L P/C 7343015**

### Proclean

3



#### Washer/Decontaminator Detergent

Proclean is a concentrated soft-metal safe detergent formulated to handle heavy duty cleaning in all water hardness conditions. Proclean effectively removes heavy soils from all types of wares without leaving a residue and complies with AS/NZS 4187:2003 and AS/NZS 4185:2001.

Proclean is available in an encapsulated, iconed and numbered bottle.

**PACK SIZE: 2X5L P/C 7683205**

### Access

3



#### Washer/Decontaminator Detergent

Access has been designed to remove heavy soils from glass, plastic and metallic instruments and is suitable for washers and decontaminators with proprietary dispensing systems. Access also complies with AS/NZS 4187:2003 and AS/NZS 4185:2001. Access is highly concentrated and gives excellent results.

Access is available in a wide-neck, colour-coded, iconed and numbered bottle.

**PACK SIZE: 2X5L P/C 7738205**

**PACK SIZE: 15L P/C 7738015**

## Automatic Rinsing

### All Dry

2



#### Washer/Decontaminator Rinse Aid

All Dry is a super-concentrated rinse additive designed for fast throughput instrument washers and decontaminators. Glass, plastic and metal surgical instruments and equipment dry quickly and spot-free. All Dry is effective in all water hardness conditions. All Dry complies with AS/NZS 4187:2003 and AS/NZS 4185:2001.

All Dry is available in an encapsulated colour-coded, iconed and numbered bottle.

**PACK SIZE: 2X5L P/C 7689205**

### Free Dry

2



#### Free Dry Super Concentrated Rinse Additive

Free Dry is a super-concentrated rinse additive designed for fast throughput instrument washers and decontaminators. Glass, plastic and metal surgical instruments and equipment dry quickly and spot-free. Free Dry is effective in all water hardness conditions. Free Dry is suitable for washers and decontaminators with proprietary dispensing systems and complies with AS/NZS 4187:2003 and AS/NZS 4185:2001. Free Dry is the same formulation as All Dry, available in a wide-neck, colour-coded, iconed and numbered bottle.

**PACK SIZE: 2X5L P/C: 7208205**

**PACK SIZE: 15L P/C: 7208015**

### Duo Rinse

2



#### Iron and Hardness Inhibiting Rinse Aid

DuoRinse was designed to prevent brown or white spotting that occurs in areas where there are problematic levels of water minerals and iron. DuoRinse is effective in all water hardness conditions. DuoRinse is designed for use in washers running within Australian Standards AS/NZS 4187:2003 and AS/NZS 4185:2001.

DuoRinse is available in a wide-neck or encapsulated iconed and numbered bottle.

**PACK SIZE: 2X5L ENCAPSULATED P/C 7347205**

**PACK SIZE: 15L P/C 7349015**

**PACK SIZE: 2X5L WIDE-NECK P/C 7349205**

## Remedial Cleaning

### Final

11



#### Remedial Acid Cleaner

Final is a remedial acidic cleaner designed to remove water scale, rust and stains from stainless steel instruments, instrument washers, stainless steel carts and trays, and sterilisers.

Final complies with AS/NZS 4187:2003 and AS/NZS 4185:2001.

**PACK SIZE: 2X5L P/C 7230205**