



Waterbath accessories

Aqua-clear water conditioner

Protects water baths and other water systems against the growth of algae, bacteria, molds and biofilms.

Just 2 mL per liter of water adds a blue tint which will fade after about 3 months reminding you to change the water and add new conditioner

- Protects pumps and thermostats and helps prevent tubing from clogging
- Permits longer use of water used in developer for color or X-ray films; will not harm film
- Can be added to the water inserted in the Water Jackets of CO₂ and regular incubators.
- Does not contain any substance listed by California's safe drinking water and toxic Enforcement Act of 1986 (Proposition 65)



Item Code	Description	Units	Price
BAT17093.0000	Cleanware aqua-clear water conditioner	100 mL	\$155.00

Lead ring flask weight

Reduces breakage, floating and spills.

- Stabilise glassware immersed in water baths; reduces floating, breakage and chemical spills
- Bright red Vikem[®] vinyl coating
- Chemical-resistant
- Will not scratch flasks

Item Code	Description	Units	Price
BAT18307.0003	Lead ring vikem covered fits flask 5-10 mL	EA	\$49.50
BAT18307.0005	Lead ring vikem covered fits flask 125-500 mL	EA	\$40.70
BAT18307.0010	Lead ring vikem covered fits flask 250-1000 mL	EA	\$70.40
BAT18307.0015	Lead ring vikem covered fits flask 500-2000 mL	EA	\$81.40
BAT18307.0020	Lead ring vikem covered fits flask 1000-4000 mL	EA	\$88.00



Microcentrifuge floating rack

Incubate up to twenty-four 1.5 mL microcentrifuge tubes in a water bath.

- Easy to load
- Microtubes are held snugly and then released when the rack is pressed down on a benchtop
- Dimensions: 11.9 x 7.8 cm



Item Code	Description	Units	Price
BAT18876.0000	Microcentrifuge floating rack 24 places for 1.5 mL tubes	Case of 4	\$82.75

Waterproof thermometers

Water-resistant/shockproof case provides flawless readings in freezers, water baths, incubators, tanks.

- Stainless-steel probe efficiently monitors temperatures in liquids, semi-solids, air/gases
- One-key operation quickly recalls highest/lowest readings over any time period
- Stainless steel probe on 120" cable with suction cap and velcro accuracy $\pm 1^{\circ}\text{C}$ between -20 and 100°C



Item Code	Description	Units	Price
FSH15-077-9E	Waterproof thermometer -50 to 300°C	EA	\$75.00

General purpose water baths

- Small foot print
- Advanced microprocessor controller, ability to store four presets
- Drain hose connection located at front of baths for easy connection
- Baths supplied with clear polycarbonate lid diffuser tray and drain hose baths 10 litre and above.
- Available in shaking water baths and recirculating water baths

Item code	Capacity	Temperature Range	Temperature Stability/ Uniformity @70°C	Work Area (L x W x D)	Overall dimensions (L x W x D)
BARTSGP02	2 litre	Ambient to 100°C	+/-0.1°C +/-0.2°C	138 x 155 x 150mm	230 x 199 x 233mm
BARTSGP05	5 litre	Ambient to 100°C	+/-0.1°C +/-0.2°C	154 x 300 x 150mm	246 x 355 x 232mm
BARTSGP10	10 litre	Ambient to 100°C	+/-0.1°C +/-0.2°C	301 x 330 x 150mm	393 x 383 x 233mm
BARTSGP20	20 litre	Ambient to 100°C	+/-0.1°C +/-0.2°C	297 x 500 x 200mm	392 x 555 x 233mm
BARTSGP28	28 litre	Ambient to 100°C	+/-0.1°C +/-0.2°C	297 x 500 x 200mm	392 x 555 x 282mm



For pricing and more information contact nzinfo@thermofisher.com

Order Placement:

For customer service, call 0800-933-966
To fax an order, use 0800-329-246
To email an order, ordersnz@thermofisher.com

Contact Us:

For customer service, call 0800-933-966
For service and calibration, call 0800-933-966
To order online: thermofisher.co.nz

Find out more at thermofisher.co.nz

ThermoFisher
SCIENTIFIC