



# March into Autumn

with great deals

**ThermoFisher**  
SCIENTIFIC

# Thermo Scientific™ Nunc™ cell culture consumables

- 10<sup>-6</sup> sterility assurance level (SAL) as per ISO 11137
- Each batch is tested with 4 cell types
  - 2 cell lines and 1 primary cell culture for monolayer formation
  - 1 cell line for cloning efficiency
- Certified non-pyrogenic
- EasYFlask™ has ergonomically designed 1/3 turn cap
- Multiwell plates have raised well rims to lower risk of cross contamination



## Quality from start to finish

Reliability • Reproducibility • Scalability • Traceability • Security • Simplicity

Item Code	Description	Unit	1 case	3+ cases
<b>Cell culture multiwell plates - Nunclon Delta™ cell culture surface - individually wrapped sterile</b>				
NUN140675	Multidish 6 well flat bottom, clear with lid	Case of 75	\$105.00	\$95.00
NUN150628	Multidish 12 well flat bottom, clear with lid	Case of 75	\$109.95	\$99.00
NUN142475	Multidish 24 well flat bottom, clear with lid	Case of 75	\$105.00	\$95.00
NUN150687	Multidish 48 well flat bottom, clear with lid	Case of 75	\$125.00	\$115.00
NUN167008	MicroWell™ 96 well flat bottom, clear with lid	Case of 50	\$58.00	\$52.00
<b>Cell culture flasks - Nunclon Delta cell culture surface - sterile</b>				
NUN156367	25cm <sup>2</sup> EasYFlask with filter cap	Case of 200	\$198.50	\$184.00
NUN156499	75cm <sup>2</sup> EasYFlask with filter cap	Case of 100	\$155.00	\$138.00
NUN159910	175cm <sup>2</sup> EasYFlask with filter cap	Case of 30	\$85.00	\$73.00

## Finnpipette™ F1 good laboratory pipetting (GLP) kits

Everything you want in your everyday, go-to pipette

- Silent, effortless volume adjustment
- Lightweight design with lowered pipetting forces reduces risk of RSI
- 5 year warranty



Item Code	Description	Unit	Price
FNP4700850N	Kit 1: 1-10 µL micro, 10-100 µL and 100-1000 µL pipettes	EA	\$790.00
FNP4700860N	Kit 2: 0.2-2 µL micro, 2-20 µL, 20-200 µL and 100-1000 µL pipettes	EA	\$890.00
FNP4700865N	Kit 3: 10-100 µL, 100-1000 µL and 1-10 mL pipettes	EA	\$790.00
FNP4701060N	Kit 4: 2-20 µL, 20-200 µL and 100-1000 µL pipettes	EA	\$790.00

# 15% off list price on all stepper tips

## Finnpipette Stepper Pipette

- Easy to use stepper pipette up to 45 dispensing
- Light weight ergonomic design
- Simple easy to set volume adjustment
- Volume range 10-5000 µL

### Stepper Pipette

Item Code	Description	Unit	Price
FNP4540000	Finnpipette stepper 10-5000 µL	EA	\$420.00
<b>Stepper tips - sterile</b>			
FNP9404173	Finntip for Finnpette stepper 0.5 mL sterile	50pcs	
FNP9404183	Finntip for Finnpette stepper 1.25 mL sterile	50pcs	
FNP9404193	Finntip for Finnpette stepper 2.5 mL sterile	50pcs	
FNP9404203	Finntip for Finnpette stepper 5 mL sterile	25pcs	
FNP9404213	Finntip for Finnpette stepper 12.5 mL sterile	25pcs	
FNP9404223	Finntip for Finnpette stepper 25 mL sterile	10pcs	
FNP9404233	Finntip for Finnpette stepper 50 mL sterile	10pcs	
<b>Stepper tips - non-sterile</b>			
FNP9404170	Finntip for Finnpette stepper 0.5 mL	100pcs	
FNP9404180	Finntip for Finnpette stepper 1.25 mL	100pcs	
FNP9404190	Finntip for Finnpette stepper 2.5 mL	100pcs	
FNP9404200	Finntip for Finnpette stepper 5 mL	50pcs	
FNP9404210	Finntip for Finnpette stepper 12.5 mL	50pcs	
FNP9404220	Finntip for Finnpette stepper 25 mL	20pcs	
FNP9404230	Finntip for Finnpette stepper 50 mL	10pcs	

\*Contact customer service on 0800-933-966 for Finntip prices



## Thermo Scientific™ Heratherm™ compact microbiological incubators

Ideal for a personalized workspace or crowded labs with small sample volumes

- Temperatures at or below ambient: 17° to 40°C with high temperature accuracy
- Save space with extremely small footprint
- 2 year warranty



Item Code	Description	Unit	Price
HEA50125882	IMC18 Heratherm compact incubator	EA	\$1,350.00

# Dependable reproducible results, flexible, easy use

Thermo Scientific mixing and temperature control drybaths

## Maximise lab space and versatility with our compact programmable drybaths

- Choice of heating/mixing, heating/cooling or heating/mixing cooling
- Excellent temperature uniformity across the block
- Easily interchangeable blocks for tubes and plates
- Digital control for easy setup and current status of run
- 2 year warranty



Item Code	Description	Unit	Price
BAR88880030	Heating cooling drybath - blocks additional	EA	\$2,595.30
BAR88880028	Heating shaking drybath - blocks additional	EA	\$2,623.40
GRA13687717	Heating shaking cooling dry bath with 1.5 ml block	EA	\$3,478.10
GRA13687721	Heating shaking cooling dry bath with 96-well microplate block	EA	\$3,478.10

## LabServ® polypropylene plasticware with moulded or printed graduation



Item Code	Description	Unit	Price
<b>Beakers</b>			
LBS80301-03	Beaker 250 mL PP low form moulded graduations	EA	\$6.50
LBS80301-05	Beaker 1000 mL PP low form moulded graduations	EA	\$12.00
LBS80301-06	Beaker 2000 mL PP low form moulded graduations	EA	\$28.50
<b>Measuring cylinders</b>			
LBS80303H-05	Measuring cylinder 250 mL PP moulded graduations	EA	\$9.00
LBS80303H-06	Measuring cylinder 500 mL PP moulded graduations	EA	\$12.50
<b>Jugs</b>			
LBSJUG3000	Jug 3 L PP printed graduations autoclavable	EA	\$25.00
LBSJUG5000	Jug 5 L PP printed graduations autoclavable	EA	\$32.50



KONICA MINOLTA

## SPAD 502 Plus chlorophyll meter

A lightweight and water resistant hand-held meter for lab and field, measuring the chlorophyll content of leaves without causing damage to the plant.

Measurements are taken by just closing the measuring head on the leaf. Since the leaf is not cut off or damaged in any other way, the same leaf can be measured as the plant grows.



Item Code	Description	Unit	Price
MINA1RT-206	SPAD-502 Plus chlorophyll meter	EA	\$3,289.00

## Spectrophotometer CM-5

Stand alone full-featured spectrophotometer with maximum flexibility and versatility that covers all your applications

15% OFF LIST PRICE

### Features and Benefits:

- Extremely easy operation with automatic calibration
- Top-port concept for Reflectance measurements as well as large transmission chamber for liquids and solids
- “Operation Wizard” with real time on screen user guidance and storage of user settings and data on USB stick



## Compact Chroma meters / colorimeters

Konica Minolta’s color difference readers are portable, battery powered, and hand-held instrument for quick, precise color control.

Item Code	Description	Unit	Price
MINA80F-107	Colour reader CR-10 Plus	EA	\$4,650.00
MINA80F-217	Colour reader CR-20	EA	\$5,348.00

\*Software included for CR-10 Plus and CR-20 readers



# Stock up on lab bottles

## Fisherbrand™ graduated laboratory bottles

- Polypropylene screw caps and pouring ring
- Borosilicate glass
- Autoclavable
- GL 45 screw thread
- Permanent white enamel marking spot



Item Code	Description	Unit	Price
FSBFB33144	Laboratory bottle 100 mL	EA	\$3.75
FSBFB33145	Laboratory bottle 250 mL	EA	\$3.95
FSBFB33146	Laboratory bottle 500 mL	EA	\$5.45
FSBFB33147	Laboratory bottle 1 L	EA	\$7.75
FSBFB33148	Laboratory bottle 2 L	EA	\$15.75
FSBFB33180	Laboratory bottle 5 L	EA	\$75.00
FSBFB33181	Laboratory bottle 10 L	EA	\$225.00
FSBFB33182	Laboratory bottle 20 L	EA	\$475.00

# Time to clean up

Clean up, clear out, declutter and rejuvenate your bench and lab space. Let us help you out with our best cleaners for the job!



Item Code	Description	Application	Unit	Price
DEC90/5	Decon 90 liquid concentrate	Manual cleaning	5 L	\$150.00
DECNEUTRA/5	Neutracon liquid	Manual cleaning	5 L	\$150.00
RBS35CONC	RBS35 concentrated cleaning agent		5 L	\$105.00
BSPSTE70/R	Ethanol denatured alcohol 70% sterile spray	General disinfection	6 x 500 mL	\$80.00
DIVHH12291	Pyroneg liquid cleaner	Manual	2 L	\$85.00
DIVHH13231	Pyroneg powder	Automatic cleaning	3 kg	\$120.00
VRKVIRKON5G	Virkon S	Virucidal disinfection	50 x 5 g tablets	\$100.00
VRKVIRKON5	Virkon S broad spectrum	Virucidal disinfection	500 g	\$70.00

# High quality Great value

## Fisher Scientific™ accumet™ pH meters



### Fisher Scientific accumet AE-150 pH Meter with Stand

Budget friendly, accurate and reliable basic pH meter

- Available as meter and stand only or in 3 different kits that include electrodes and temperature probes
- BNC input for electrode compatibility
- Automatic or manual temperature compensation
- pH slope and offset display
- One to five point pH calibration and mV offset adjustment for ORP
- 3 Years Warranty

Item Code	Description	Unit	Price
ACC13636-AE150	AE150 pH meter with stand	EA	\$596.10
ACC13636-AE151	AE150 pH meter - education set	EA	\$682.90
ACC13636-AE152	AE150 pH meter - bio set	EA	\$688.20
ACC13636-AE153	AE150 pH meter - 3 in one set	EA	\$699.30



### Fisher Scientific accumet AB series meters

All models have the below features and are available as meter only or in kits that include electrodes and temperature probes

- AB-150: Able to measure pH/mV/Temperature
- AB-200: Able to measure pH/mV/Temperature/Conductivity/TDS/Resistivity/Salinity
- AB-250: Able to measure pH/mV/Temperature/Ion Selective Measurements
- 0.001 pH resolution and  $\pm 0.002$  pH accuracy
- Date and time logging to meet GLP requirements
- BNC input for electrode compatibility and RS-232 output for "21 CFR Part 11 Compliance"
- Password protection and 500 point memory with manual or automatic interval logging
- 3 Years Warranty

Item Code	Description	Unit	Price
ACC13636-AB150A	accumet AB150 pH meter only	EA	\$918.70
ACC13636-AB150	accumet AB150 pH kit	EA	\$1,018.00
ACC13636-AB150B	accumet AB150 pH bio kit	EA	\$1,088.80
ACC13636-AB200A	AB200 pH/conductivity meter only	EA	\$1,385.20
ACC13636-AB200	accumet AB200 pH/conductivity kit	EA	\$1,593.50
ACC13636-AB250A	accumet AB250 pH/ISE meter	EA	\$1,058.80
ACC13636-AB250	accumet AB250 pH/ISE kit	EA	\$1,118.00



Don't forget

We have a full range of chemicals,  
check it out online at  
[www.thermofisher.co.nz/chemicals](http://www.thermofisher.co.nz/chemicals)

Don't forget

to check out more deals in our  
lab essential promotion  
[www.thermofisher.co.nz/labessentials](http://www.thermofisher.co.nz/labessentials)

---

**In New Zealand:**

For customer service, call 0800-933-966

To fax an order, use 0800-329-246

To email an order, [ordersnz@thermofisher.com](mailto:ordersnz@thermofisher.com)

Find out more at:

[thermofisher.co.nz/promotions](http://thermofisher.co.nz/promotions)

**ThermoFisher**  
SCIENTIFIC

Offers valid till 31st March 2018. All offers subject to Thermo Fisher standard terms and conditions. Errors and omission exempted.  
Specifications, terms and pricing are subject to change. Please consult your local sales representative for details. © 2018 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. 1519099635