



Kimble commonly used glassware collection

Highest quality borosilicate glassware for your lab

Beakers

Low Form

KIMAX® beakers offer excellent mechanical strength and durability, while providing high resistance to chemical attack and thermal shock. They have been a staple in research laboratories for many generations.

Product Code	Description	Price	Product Code	Description	Price
KIM14000-10	Beaker Glass 10 mL	\$6.20	KIM14000-400	Beaker Glass 400 mL	\$7.00
KIM14000-20	Beaker Glass 20 mL	\$4.20	KIM14000-600	Beaker Glass 600 mL	\$8.10
KIM14000-30	Beaker Glass 30 mL	\$5.20	KIM14000-800	Beaker Glass 800 mL	\$10.75
KIM14000-50	Beaker Glass 50 mL	\$5.20	KIM14000-1000	Beaker Glass 1000 mL	\$15.25
KIM14000-100	Beaker Glass 100 mL	\$5.60	KIM14000-1500	Beaker Glass 1500 mL 16pcs	\$286.00
KIM14000-150	Beaker Glass 150 mL	\$6.00	KIM14000-2000	Beaker Glass 2000 mL	\$29.60
KIM14000-250	Beaker Glass 250 mL	\$6.70	KIM14000-4000	Beaker Glass 4000 mL 6pcs	\$540.00



Low Form Heavy Duty

KIMAX® heavy duty beakers offer superior mechanical strength and durability. Improved safety when used under extreme conditions such as mechanized washing and rough handling.

Product Code	Description	Price
KIM14005-250	Beaker Glass 250 mL	\$12.30
KIM14005-400	Beaker Glass 400 mL	\$16.00
KIM14005-600	Beaker Glass 600 mL	\$18.70
KIM14005-1000	Beaker Glass 1000 mL	\$35.50
KIM14005-2000	Beaker Glass 2000 mL	\$72.50
KIM14005-4000	Beaker Glass 4000 mL	\$135.00



Tall Form

KIMAX® Berzelius beakers offer excellent mechanical strength and durability, while providing high resistance to chemical attack and thermal shock. Ideal for use when performing titrations.

Product Code	Description	Price
KIM14030-100	Beaker Glass 100 mL spouted	\$7.95
KIM14030-200	Beaker Glass 200 mL spouted	\$8.05
KIM14030-300	Beaker Glass 300 mL spouted	\$9.00
KIM14030-400	Beaker Glass 400 mL spouted	\$10.50
KIM14030-600	Beaker Glass 600 mL spouted	\$14.00
KIM14030-1000	Beaker Glass 1000 mL spouted	\$18.50



All glassware is borosilicate 3.3, conforming to ASTM E348, Type 1, Class A

Our glassware catalogue covers a range of essential glassware for every lab.



View our catalogue online or contact us on 0800-933-966 for your hard copy.

thermofisher.co.nz/glassware



For added convenience, the Kimble Chase catalogue is also available online through the link below.

thermofisher.co.nz/kimblechase

Media bottles

GL45 Media Bottles

Ideal for general laboratory use including mixing, storing or transporting culture media, chemicals or solvents. Borosilicate 3.3 glass, graduated, with PP cap and pouring ring.

Product Code	Description	Price
KIM14395-100	Bottle lab 100 mL	\$13.00
KIM14395-250	Bottle lab 250 mL 10 pcs	\$144.00
KIM14395-500	Bottle lab 500 mL 10 pcs	\$175.00
KIM14395-1000	Bottle lab 1000 mL	\$21.00
KIM14395-2000	Bottle lab 2000 mL	\$55.00
KIM14395-5000	Bottle lab 5000 mL	\$181.00
KIM14395-10000	Bottle lab 10000 mL	\$315.00



*All bottles come with PP screw cap and pouring ring and are autoclavable. Borosilicate glass.

Product Code	Description	Price
KIM1435P-45	Pouring ring clear GL45 PP autoclavable 10pcs	\$22.00
KIM14395E-452	Pouring ring red GL45 ETFE autoclavable 10pcs	\$45.00
KIM14395C-45	Screw cap blue GL45 PP autoclavable 10pcs	\$60.00
KIM14395H-452	Screw cap red GL45 PBT with ETFE protected seal 10pcs	\$220.00




Raysorb®

Designed to protect contents from UV rays; ideal for light-sensitive products. Amber graduated borosilicate bottles with PP cap and pouring ring.

Product Code	Description	Price
KIM14399-100	Bottle lab 100 mL	\$60.00
KIM14399-250	Bottle lab 250 mL	\$92.80
KIM14399-500	Bottle lab 500 mL	\$87.80
KIM14399-1000/EA	Bottle lab 1000 mL	\$96.00
KIM14399-2000	Bottle lab 2000 mL	\$142.00
KIM14399-5000	Bottle lab 5000 mL	\$473.00
KIM14399-10000	Bottle lab 10000 mL	\$705.00




Raysorb: Better than Amber



Kimcote® (Plastic Coated)

Product Code	Description	Price
KIMKC14395-100	Bottle lab 100 mL 4pcs	\$400.00
KIMKC14395-250	Bottle lab 250 mL 4pcs	\$460.00
KIMKC14395-500	Bottle lab 500 mL 4pcs	\$540.00
KIMKC14395-1000	Bottle lab 1000 mL 4pcs	\$800.00
KIMKC14395-2000	Bottle lab 2000 mL 4pcs	\$1,130.00
KIMKC14395-5000	Bottle lab 5000 mL	\$795.00
KIMKC14395-10000	Bottle lab 10000 mL	\$1,220.00



Erlenmeyer Flasks

Narrow Neck

KIMAX® borosilicate flasks for economy and versatility. These flasks are the choice for general laboratory usage. ASTM quality heavy tooled necks for strength.



Product Code	Description	Price
KIM26500-25	Flask erlenmeyer 25 mL 48pcs	\$370.00
KIM26500-50	Flask erlenmeyer 50 mL	\$7.20
KIM26500-125	Flask erlenmeyer 125 mL 48pcs	\$308.00
KIM26500-250	Flask erlenmeyer 250 mL	\$7.70
KIM26500-300	Flask erlenmeyer 300 mL 48pcs	\$458.00
KIM26500-500	Flask erlenmeyer 500 mL	\$11.80
KIM26500-1000	Flask erlenmeyer 1000 mL 24pcs	\$375.00
KIM26500-2000	Flask erlenmeyer 2000 mL	\$36.60
KIM26500-4000	Flask erlenmeyer 4000 mL	\$115.80
KIM26500-6000	Flask erlenmeyer 6000 mL	\$210.00

Wide Neck

Product Code	Description	Price
KIM26650-125	Flask erlenmeyer 125 mL	\$8.60
KIM26650-250	Flask erlenmeyer 250 mL	\$9.20
KIM26650-500	Flask erlenmeyer 500 mL 36pcs	\$416.00
KIM26650-1000	Flask erlenmeyer 1000 mL 24pcs	\$462.00
KIM26650-2000	Flask erlenmeyer 2000 mL 8pcs	\$288.00



Buy online with goDirect to reduce cost & improve efficiency

- Stock and Pricing for your account.
- Save templates, save time.
- See your backorders and ETAs.



Order Placement:

For customer service, call 0800-933-966
To fax an order, use 0800-329-246
To email an order, ordersnz@thermofisher.com

Contact Us:

For customer service, call 0800-933-966
For service and calibration, call 0800-933-966
To order online: thermofisher.co.nz/godirect

Shop online at thermofisher.co.nz/godirect

ThermoFisher
SCIENTIFIC

Terms & conditions: Offers valid till 30th September 2018. This promotion is open only to customers in Australia and New Zealand. All offers subject to Thermo Fisher standard terms and conditions. Errors and omission exempted. Specifications, terms and pricing are subject to change. Please consult your local sales representative for details. © 2018 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. Kimble is a trademark of DWK Life Sciences. **1520550493**