



ART Barrier Tip Reload System

The green alternative to barrier tip packaging

The Thermo Scientific™ ART™ barrier tips offer a complete selection of tips to fit virtually any pipette or application in your laboratory. Available in the versatile hinged and reload systems.

ART Hinged and Reload Rack Benefits

- Hinged rack offers one handed access to tips – no more tips spilling on the bench or floor maximising productivity.
- Individually packaged reload blister tubs easily insert into the hinged rack.
- ART reload products uses **60% less plastic** than a racked product.
- All packaging plastic is completely recyclable.
- **Save 44%** storage space using the reload system.

ART Reload System Steps



Step 1

Hold the top of the blister tub and peel the paper on the bottom of the tub.



Step 2

Align the blister tub with the hinged rack.



Step 3

Push down directly, locking the blister tub into the hinged rack.

Environmentally friendly, Space saving design, Hinged covers and Efficient reloading

ART Barrier Tips Features

- ART tips deliver maximum sample recovery through the use of low-retention technology.
- The self-sealing barrier provides 100% protection against contamination from aerosols and liquids.
- All tips are certified as RNase, DNase and pyrogen free.
- Available in the versatile hinged rack and reload system

Special Offer

Buy a pack of reload inserts and receive a pack of hinged rack tips for half price!

Use promo-code: **ARTPROMO**

Item Code -Reload Insert	Description	Item Code – Hinge Rack
MBP2139-05-RI	Barrier tips, low retention, reload insert, 10 μ L, 96 tips x 10 Racks	MBP2139-05-HR
MBP2140-05-RI	Barrier tips, low retention, reload insert, extended length, 10 μ L, 96 tips x 10 Racks	MBP2140-05-HR
MBP2149E-05-RI	Barrier tips, low retention, reload insert, 10 μ L, 96 tips x 10 Racks	MBP2149-05-HR
MBP2065E-05-RI	Barrier tips, low retention, reload insert, 100 μ L, 96 tips x 10 Racks	MBP2065E-05-HR
MBP2069-05-RI	Barrier tips, low retention, reload insert, 200 μ L, 96 tips x 10 Racks	MBP2069-05-HR
MBP2179-05-RI	Barrier tips, low retention, reload insert, 1000 μ L, 96 tips x 8 Racks	MBP2179-05-HR

Order Placement:

For customer service, call 0800-933-966
To fax an order, use 0800-329-246
To email an order, ordersnz@thermofisher.com

Contact Us:

For customer service, call 0800-933-966
For service and calibration, call 0800-933-966
To order online: thermofisher.co.nz

Find out more at thermofisher.co.nz