



Vintage Essentials Catalogue

Your source for essential instrumentation
and consumables for wine

ThermoFisher
SCIENTIFIC

Essential lab supplies for this year's vintage

From the largest wine labs to the simple pH strips or titration set up, Thermo Fisher Scientific has you covered for general and specialty labware. We have over 100,000 products listed on our online store, with stock available from hundreds of the best agencies for wine analysis.

We have endeavoured to include the most relevant products in this publication. From equipment to consumables to enzyme kits, we have the largest portfolio of products and brands to offer you.

Browse here or online, contact your account manager or customer service, or place an order via our national customer service centre on **1300 735 292**, alternatively visit **thermofisher.com.au/godirect**

We look forward to helping stock your lab with the best laboratory range available.

brands we supply

thermo
scientific

invitrogen

applied
biosystems



ACRÖS
ORGANICS



KIMBLE
CHASE



PYREX®

LabServ®

sartorius

BS Bellingham
+ Stanley
a xylem brand

Fisherbrand

Whatman®
Part of GE Healthcare

ELGA

MERCK

Tuttnauer
Your Sterilization & Infection Control Partners

Contents



Chemicals & Cleaners

4



Consumables & General Labware

10



Equipment

28



Filtration

42



Microbiology & Hygiene Monitoring

46



Reagent Kits & Analysers

50



Cellular & Molecular Tools for Wine Analysis

60



Chemicals & Cleaners

Acids	5
Dries	5
Solutions	6
Solvents	7
Standards	8
Buffers	9
Test Strips	9
Cleaners	9

Full range including dry chemicals, solutions, solvents, acids, volumetric solutions, biochemicals, standards, high purity solvents & acids and specialty chemicals.

Acids

From the competitively priced LabServ Pronalys and Ajax Finechem acids to the superior analytical grade Fisher Chemical acids, you are bound to be pleased with the quality of the chemicals below - commonly used in the wine industry.

Item	Description	Unit
AJA1-2.5LPL	Acetic acid glacial UNIVAR	2.5L
FSBF/1900/PB17	Formic acid CertiFied AR 98-100%	2.5L
BSPVL701.5	Hydrochloric acid 0.1M (0.1N) Volusol plastic carboy	5L
BSPVL700.5	Hydrochloric acid 1.0M (1.0N) Volusol plastic carboy	5L
AJA256-2.5LPL	Hydrochloric acid 32% UNIVAR	2.5L
AJA1367-2.5LGL	Hydrochloric acid 36% UNIVAR	2.5L
AJA341-2.5LGL	Nitric acid 70% UNIVAR ACS	2.5L
AJA371-2.5LPL	Orthophosphoric acid 85% UNIVAR ACS	2.5L
AJA1530-2.5LPL	Orthophosphoric acid 25% UNIVAR	2.5L
AJA534-2.5LGL	Sulphuric acid 98% UNIVAR ACS	2.5L



Dries

Thermo Fisher Scientific have a wide range of chemical dries and salts including powders and granules. These are available in packs of 5g-5kg.

Item	Description	Unit
BSPAL847.250	Anti bumping granules	250g
BSPBL952.500	Boric acid AR ACS	500g
ACR151320050	Bromocresol green sodium salt water soluble	5g
FSBC/2360/53	Calcium sulphate dihydrate CertiFied AR	500g
AJA125-500G	Calcium carbonate UNIVAR	500g
AJA127-500G	Calcium chloride dihydrate UNIVAR	500g
AJA977-500G	Charcoal activated decolourising TECHNICAL	500g
BSPCL942.500	Copper (II) sulphate pentahydrate Pronalys AR ACS	500g
AJA478-500G	Di-sodium hydrogen orthophosphate AR	500g
BSPSL945.500	EDTA disodium salt dihydrate Pronalys AR ACS	500g
FSBJ/3700C/05	EDTA trisodium salt 0.1M Solutrate concentrate	1 ampoule
FSBA/8882/53	L-Ascorbic acid CertiFied AR	500g
FSBI/0505/50	Iodine indicator for wine titrations SpeciFied LR	250g
FSBM/4870/46	Methylene blue C.I. 52015 redox indicator	25g



GL= Glass PL= Plastic

Dries (Continued)

Item	Description	Unit
ACR160280250	N-Acetyl-L-cysteine 98%	25g
FSBP/3298/44	o-Phthalaldehyde laboratory Reagent grade	5g
BSPPT403.100	Phenolphthalein powder R Grade pH indicator	100g
FSBP/3700/48	Polyvinyl pyrrolidone (P.V.P.) GP	100g
BSPPL402.500	Potassium chloride Pronalys AR ACS	500g
AJA409-500G	Potassium iodide UNIVAR	500g
BSPPS303.500	Potassium sodium (+)-tartrate tetrahydrate Pronalys AR	500g
BSPSL519.1	Silica gel 2-6mm self indicating orange	1kg
AJA8745-500G	Silica gel 2-5mm self indicating orange	500g
BSPSL969.500	Sodium carbonate anhydrous Pronalys AR ACS	500g
FSBS/2920/53	Sodium carbonate anhydrous CertiFied AR	500g
BSPSL944.5	Sodium chloride Pronalys AR ACS	5kg
BSPSL944.500	Sodium chloride Pronalys AR ACS	500g
BSPSL900.500	Sodium hydrogen carbonate Pronalys AR ACS	500g
BSPSL740.5	Sodium hydroxide pellets Pronalys AR ACS	5kg
BSPSL740.500	Sodium hydroxide pellets Pronalys AR ACS	500g
FSBS/6640/53	Sodium sulphate anhydrous CertiFied AR	500g
FSBT/P630/60	Tris (hydroxymethyl) methylamine 99.8+%	1kg
FSBBP152-500	Tris Base Molecular Biology grade	500g
BSPSL808.500	Sodium thiosulphate Pronalys AR ACS	500g
BSPST730.500	Starch soluble Pronalys AR	500g
FSBS/8600/60	Sucrose CertiFied AR	1kg



Solutions

Available in a broad range of volumes, purities and concentrations, these Fisher Chemical, LabServ Pronalys and Ajax Finechem solutions are ideal for use in your winery.

Item	Description	Unit
FSBA/0080/08	Acetaldehyde SpeciFied LR	500mL
BSPAL724.2.5	Acetone Pronalys AR ACS glass bottle	2.5L
AJA43-2.5LGL	Ammonia solution 28 - 30% UNIVAR	2.5L
BSPWL012.5	Distilled water	5L
BSPWL012.20	Distilled water	20L
FSBB/4320L/08	Bromocresol green 0.04% solution pH 3.6-5.2 yellow to blue	500mL
AJA799-2.5LPL	Fehlings solution No.1 LABCHEM	2.5L
AJA800-2.5LPL	Fehlings solution No.2 LABCHEM	2.5L
BSPA5.2.5	Hydrogen peroxide 30% 100/vol Pronalys AR ACS	2.5L
BSPA5.500	Hydrogen peroxide 30% 100/vol Pronalys AR ACS	500mL
AJA1868-100ML	Rankin indicator solution LABCHEM	100mL



Solutions (Continued)

Item	Description	Unit
FSBH/1800/15	Hydrogen peroxide >30% w/v CertiFied AR	1L
ACR124220010	Iodine 0.05M 0.1N Volumetric solution	1L
BSPPL460.100	Phenolphthalein 1% solution	100mL
BSPVL703.250	Potassium chloride 3M Volusol plastic bottle	250mL
BSPVL711.2.5	Sodium hydroxide 0.01N Volusol plastic bottle	2.5L
BSPVL713.2.5	Sodium hydroxide 0.1N Volusol plastic bottle	2.5L
BSPVL712.2.5	Sodium hydroxide 1.0N Volusol plastic bottle	2.5L
AJA2539-2.5LPL	Sodium hydroxide 0.333M UNIVOL	2.5L
AJA485-5L	Sodium hypochlorite solution 12.5% TECHNICAL	5L
FSBJ/7950C/90	Sodium thiosulphate 0.1mol/L conc. makes 1L	6pcs
ROA2653/2.5L	Z1 solution Rebelein sugar	2.5L
ROA2654/2.5L	Z2 solution Rebelein sugar	2.5L
ROA2655/2.5L	Z3 solution Rebelein sugar	2.5L
ROA2656/2.5L	Z4 solution Rebelein sugar	2.5L
ROA2657/2.5L	Z5 solution Rebelein sugar	2.5L
ROA2658/2.5L	Z6 solution Rebelein sugar	2.5L
R40010	Acridine orange btl	250mL



Solvents

Use of high quality analytical grade solvents in your QC laboratory is important for ensuring reliable testing of your wines. The solvents listed are available in 500mL up to 2.5L depending on your needs, or contact us for other volumes.

Item	Description	Unit
FSBB/4850/17	Butan-1-ol CertiFied AR	2.5L
FSBD/2050/08	Diethanolamine SpeciFied extra pure	500mL
BSPEL978.2.5	Ethanol 95% v/v Pronalys AR ACS glass bottle	2.5L
BSPEL975.2.5	Ethanol 99.7-100% absolute Pronalys AR ACS glass bottle	2.5L
AJA348-2.5LGL	Iso-octane (2,2,4-Trimethylpentane) UNIVAR ACS	2.5L
AJA318-2.5LPL	Methanol UNIVAR ACS	2.5L



Standards

Item	Description	Unit
ROA0653/500ML	Aluminium standard 1,000mg/L in 2% nitric acid	500mL
ROA0654/500ML	Aluminium standard 10,000ppm in 2% nitric acid	500mL
ROA2579/500ML	Antimony standard 1,000mg/L in 2% nitric acid	500mL
ROA2610/500ML	Antimony standard 10,000mg/L in 2% nitric acid	500mL
ROA0665/500ML	Boron standard 1,000mg/L in water matrix	500mL
ROA0666/500ML	Boron standard 10,000mg/L in water matrix	500mL
ROA0669/500ML	Calcium standard 1,000mg/L in 2% nitric acid	500mL
ROA0670/500ML	Calcium standard 1,0000mg/L in 2% nitric acid	500mL
ROA0679/500ML	Copper standard 1,000mg/L in 2% nitric acid	500mL
ROA0680/500ML	Copper standard 10,000ppm in 2% nitric acid	500mL
ROA0701/500ML	Iron standard 1,000mg/L in 2% nitric acid	500mL
ROA0702/500ML	Iron standard 10,000mg/mL in 2% nitric acid	500mL
ROA0705/500ML	Lead standard 1,000mg/L in 2% nitric acid	500mL
ROA0706/500ML	Lead standard 10,000ppm in 2% nitric acid	500mL
ROA0711/500ML	Manganese standard 1,000mg/L in 2% nitric acid	500mL
ROA0712/500ML	Manganese standard 10,000ppm in 2% nitric acid	500mL
ROA0713/500ML	Magnesium standard 1,000mg/L in 2% nitric acid	500mL
ROA0714/500ML	Magnesium standard 10,000mg/L in 2% nitric acid	500mL
ROA0717/500ML	Molybdenum standard 1,000ppm in water matrix	500mL
ROA0718/500ML	Molybdenum standard 10,000mg/L in water matrix	500mL
ROA0727/500ML	Phosphorus standard 1,000mg/L in water matrix	500mL
ROA0728/500ML	Phosphorus standard 10,000mg/L in water matrix	500mL
ROA0731/500ML	Potassium standard 1,000mg/L in 2% nitric acid	500mL
ROA0732/500ML	Potassium standard 10,000mg/L in 2% nitric acid	500mL
ROA0744/500ML	Selenium standard 10,000mg/L in 2% nitric acid	500mL
ROA0747/500ML	Silver standard 1,000ppm in 2% nitric acid	500mL
ROA0749/500ML	Sodium standard 1,000mg/L in 2% nitric acid	500mL
ROA0750/500ML	Sodium standard 10,000mg/L in 2% nitric acid	500mL
ROA0764/500ML	Tin standard 10,000mg/L in 10% HCl	500mL
ROA0777/500ML	Zinc standard 1,000mg/L in 2% nitric acid	500mL
ROA0778/500ML	Zinc standard 10,000mg/L in 2% nitric acid	500mL
SPXWINE-3	Organic Wine Standard, 14 component mix 10,000ppm	1mL

Chemicals & Cleaners

Buffers

LabServ Pronalys and Fisher Chemical buffers are NIST traceable for quality assurance. Buffers are available in pH 4, 7 and 10 as well as in tablet form.

Item	Description	Unit
BSPA97.1	Buffer solution pH 4 colour-coded red	1L
BSPA97.500	Buffer solution pH 4 colour-coded red	500mL
BSPA17.500	Buffer solution pH 4 colourless	500mL
BSPA98.1	Buffer solution pH 7 colour-coded green	1L
BSPA98.500	Buffer solution pH 7 colour-coded green	500mL
BSPA18.500	Buffer solution pH 7 colourless	500mL
BSPA99.1	Buffer solution pH 10 colour-coded blue	1L
BSPA99.500	Buffer solution pH 10 colour-coded blue	500mL



Test Strips

Fisher Chemical have pH strips to cover the full pH range or targeted to the range of interest. HACH also offer water quality test strips for quick on-site evaluations of free chlorine, total chlorine, total hardness, total alkalinity and pH.

Item	Description	Unit
HAC27552-50	5-in-1 water quality test strip kit for chlorine, hardness, alkalinity and pH	50pcs
FSBFB33003	pH indicator stick non-bleeding 0-14	100pcs
FSBFB33005	pH indicator stick non-bleeding 0-6.0	100pcs
FSBFB33015	pH indicator stick non-bleeding 7.0-14	100pcs



Cleaners

Ensure your winery laboratory and consumables are kept spotless with the cleaners and detergents listed below.

Item	Description	Unit
DEC90-5L	Decon-90™ liquid concentrate. A totally efficient surface active agent for laboratory, medical & specialised industrial cleaning applications. Not suitable for use on non-ferrous metals.	5L
DIV12291	Pyronex liquid concentrate. Mildly alkaline cleaner and disinfectant for manual and ultrasonic washing.	2L
DIV13231	Pyronex powder concentrate. Mildly alkaline cleaner and disinfectant for manual and ultrasonic washing.	3kg



Item	Description	Unit
BSPSTE70/R	Ethanol denatured alcohol 70%. Sterile filtered through 0.22 micron. Trigger spray individually sealed in overwrap	6 x 500mL
BSPSTI70	Iso propyl alcohol 70%. Sterile filtered through 0.22 micron. Trigger spray individually sealed in overwrap	6 x 500mL



Consumables & General Labware

Autoclave Bags & Tape	11	Magnetic Stirrers	19
Brushes	11	Markham & Rankine Apparatus	20
Beakers	11	Measuring Cylinders	21
Bench Protectors	12	Parafilm	21
Bottles	12	Petri Dishes	21
Burettes	13	Pipettes	21
Carboys	14	Safety Glasses	24
Centrifuge Tubes	14	Slides	24
Containers	14	Spatulas	24
Coverslips	15	Syringes & Needles	24
Culture Tubes	15	Test Tube Racks	25
Cuvettes	16	Timers	25
Draining Racks	16	Thermometers	25
Flasks	16	Vials & Closures	26
Funnels	17	Wash Bottles	27
Gloves	17	Weigh Boats	27
Hydrometers	18		
Inoculating Loops & Spreaders	18		
Jugs	19		
Label Printer	19		

Our extensive capabilities in consumables and disposables from many of the world's leading brands will cover all of your laboratory needs.

Autoclave Bags & Tape

Item	Description	Unit
LBS3500/502	Autoclave bag biohazard 300 x 560mm 121°C yellow	50pcs
LBS3500/504	Autoclave bag biohazard 600 x 900mm 121°C yellow	50pcs
LBS3500/503	Autoclave bag biohazard 410 x 610mm 121°C yellow	50pcs
3M1322-12	Autoclave tape 12mm x 55m roll	EA
3M1322-18	Autoclave tape 18mm x 55m roll	EA
3M1322-24	Autoclave tape 24mm x 55m roll	EA
DRA500025	Sterilope Sterilising Satchel No. 2 paper 34x20x3.5cm	500pcs
DRA500055	Sterilope Sterilising Satchel No. 5 paper 26.5x16.5x7cm	1000pcs
DRA500067	Sterilope Sterilising Satchel No. 7 paper 27x70x3.5cm	2000pcs



Brushes

Item	Description	Unit
BEL231144	Brush for test tubes 16-25mm bristle fan end	EA
BELBT005	Brush for test tubes 20x65x250mm nylon fan tipped	EA
BEL231012	Brush bottle bristle type 2 50x75mm head 330mm length	5pcs
BELBB012	Brush bottle 70x100x550mm nylon tufted	EA
FSB03-621A	Brush nessler and cylinder bristle 10cm x 4.4cm radial double fan tip galvanised wire handle 30cm	12pcs
FSB03-614	Brush Buret bristle 8.2cm x 1.9cm wire handle galvanised to resist erosion length 91cm	12pcs



Beakers (Glass)

Heavy duty beakers and tall form beakers are also available.

Item	Description	Unit
FSBFB10050	Beaker Glass 50mL low form	EA
FSBFB100100	Beaker Glass 100mL low form	EA
FSBFB100150	Beaker Glass 150mL low form	EA
FSBFB100250	Beaker Glass 250mL low form	EA
FSBFB100400	Beaker Glass 400mL low form	EA
FSBFB100600	Beaker Glass 600mL low form	EA
FSBFB1001000	Beaker Glass 1000mL low form	EA



Consumables & General Labware

Beakers (Plastic)


Polypropylene graduated beakers. Fully autoclavable.

Item	Description	Unit
LBS80301-02	Beaker 100mL low form PP with moulded graduations	EA
LBS80301-03	Beaker 250mL low form PP with moulded graduations	EA
LBS80301-04	Beaker 500mL low form PP with moulded graduations	EA
LBS80301-05	Beaker 1L low form PP with moulded graduations	EA
LBS80301-06	Beaker 2L low form PP with moulded graduations	EA
NAL1201-0050	Beaker 50mL PP low form with printed graduations	12pcs
NAL1201-0100	Beaker 100mL PP low form with printed graduations	12pcs
NAL1201-0250	Beaker 250mL PP low form with printed graduations	6pcs
NAL1201-0600	Beaker 600mL PP low form with printed graduations	4pcs
NAL1201-1000	Beaker 1L PP low form with printed graduations	3pcs



Bench Protectors

Item	Description	Unit
NAL74018-00	Bench protector VERSI-DRY Super Absorbent Mat 18in x 20in (45.7cm x 50.8cm)	350pcs
NAL74000-00	Bench protector VERSI-DRY Super Absorbent Roll 20in x 250ft (50.8cm x 76.2m)	Roll
WHA2300-731	Benchkote reel 46cm x 50m	Roll
WHA2300-916	Benchkote sheets 46x57cm	50pcs
NAL62065-00	VERSI-DRY Standard Absorbancy roll 50.8cm x 45.7cm	2pcs



Bottles (Glass)

Fisherbrand laboratory bottles are manufactured using clear borosilicate glass with white graduations. Supplied with blue polypropylene screw cap and clear pouring ring. Fully autoclavable.

Item	Description	Unit
FSBFB33144	Bottle lab 100mL graduated with PP screw cap	EA
FSBFB33145	Bottle lab 250mL graduated with PP screw cap	EA
FSBFB33146	Bottle lab 500mL graduated with PP screw cap	EA
FSBFB33147	Bottle lab 1L graduated with PP screw cap	EA
FSBFB33148	Bottle lab 2L graduated with PP screw cap	EA
FSBFB33180	Bottle lab 5L graduated with PP screw cap	EA
FSBFB33181	Bottle lab 10L graduated with PP screw cap	EA
FSBFB33182	Bottle lab 20L graduated with PP screw cap	EA



Consumables & General Labware

Bottles (Plastic)

Item	Description	Unit
NAL2104-0001	Bottle 30mL HDPE wide mouth with PP screw cap	12pcs
NAL2104-0002	Bottle 60mL HDPE wide mouth with PP screw cap	12pcs
NAL2104-0004	Bottle 125mL HDPE wide mouth with PP screw cap	12pcs
NAL2104-0008	Bottle 250mL HDPE wide mouth with PP screw cap	12pcs
NAL2104-0016	Bottle 500mL HDPE wide mouth with PP screw cap	12pcs
NAL2104-0032	Bottle 1000mL HDPE wide mouth with PP screw cap	6pcs
LBSLBHD60/00/28N	Bottle only 60mL round HDPE natural, requires 28mm screw cap LBSLBCAP28N	600pcs
LBSLBHD125/00/28N	Bottle only 125mL round HDPE natural, requires 28mm screw cap LBSLBCAP28N	300pcs
LBSLBHD250/00/28N	Bottle only 250mL round HDPE natural, requires 28mm screw cap LBSLBCAP28N	160pcs
LBSLBHD500/00/28N	Bottle only 500mL round HDPE natural, requires 28mm screw cap LBSLBCAP28N	105pcs
LBSLBCAP28N	Cap only 28mm PP natural 28mm screw cap for LBHD bottles	4000pcs
LBSBPC150	Bottle 150mL PC natural PP screw cap	120pcs
LBSBPC300	Bottle 300mL PC natural PP screw cap	100pcs
LBSBPC600	Bottle 600mL PC natural PP screw cap	60pcs



Bottle Top Dispenser for aqueous solutions

Item	Description	Unit
GRF101.080-27	Dispenser bottle top 0.4-2mL Optifix basic	EA
GRF101.080-33	Dispenser bottle top 1-5mL Optifix basic	EA
GRF101.080-38	Dispenser bottle top 2-10mL Optifix basic	EA
GRF101.080-44	Dispenser bottle top 6-30mL Optifix basic	EA
GRF101.080-48	Dispenser bottle top 10-50mL Optifix basic	EA
GRF101.080-51	Dispenser bottle top 20-100mL Optifix basic	EA
GRF101.080-53	Dispenser bottle top 40-200mL Optifix basic	EA
GRF101.080-56	Dispenser bottle top 60-300mL Optifix basic	EA
GRF101.080-58	Dispenser bottle top 100-500mL Optifix basic	EA



Burettes (Automatic)

Automatic burette Schellbach

Item	Description	Unit
GLS118.240.01	10mL with intermediate PTFE key & needle valve stopcock	EA
GLS118.240.02	25mL with intermediate PTFE key & needle valve stopcock	EA
GLS118.240.03	50mL with intermediate PTFE key & needle valve stopcock	EA
GLS118.240.01A	10mL with needle valve stopcock no intermediate key	EA
GLS118.240.02A	25mL with needle valve stopcock no intermediate key	EA
GLS118.240.03A	50mL with needle valve stopcock no intermediate key	EA



Consumables & General Labware

Carboys

Item	Description	Unit
NAL2250-0050	Carboy 20L PP with handles and 83B PP closure, autoclavable	EA
NAL2318-0050	Carboy 20L LDPE with PP spigot (tap), handles and 83B screw cap	EA
NAL2319-0020	Carboy 10L PP with PP spigot (tap), handles and 83B screw cap, autoclavable	EA
NAL2319-0050	Carboy 20L PP with PP spigot (tap), handles and 83B screw cap, autoclavable	EA



Centrifuge Tubes (Conical)

Item	Description	Unit
NUN339650	15mL conical centrifuge tube PP sterile bulk pack	500pcs
NUN339651	15mL conical centrifuge tube PP sterile racked	500pcs
NUN339652	50mL conical centrifuge tube PP sterile bulk pack	500pcs
NUN339653	50mL conical centrifuge tube PP sterile racked	300pcs
LBSCT5203	50mL conical centrifuge tube PP non sterile white cap	500pcs
LBSCT5206	50mL centrifuge tube PP graduated separate white PE cap	500pcs



Containers (Plastic)

Produced by our LabServ plastics manufacturing plant, we have a large range of plastic containers that have numerous uses for sample collection in the field and laboratory. Please contact us for a LabServ catalogue or for more information.

Item	Description	Unit
LBS33030	Container 400mL PET natural PP screw cap square 64mm x 110mm high	100pcs
LBS33030W	Container 400mL PET natural wadded PP screw cap square	100pcs
LBS33030AW	Container 400mL PET amber natural wadded PP screw cap square	100pcs
LBS33032	Container 400mL PET label natural screw cap square	100pcs
LBS33034W	Container 400mL PET square natural wadded screw cap gamma sterile	60pcs
LBS33036W	Container 400mL PET square label wadded screw cap gamma sterile	60pcs
LBS351000	Container 1000mL PET square natural wadded screw cap PP 84mm x 171mm high	42pcs
LBS351000X	Container 1000mL PET square natural wadded screw cap PP gamma sterile	42pcs



Consumables & General Labware


Containers with Sodium Thiosulphate for Water Testing

Item	Description	Unit
LBSST120X	Container 120mL PP 10mg sodium thiosulphate spray labelled blue PE screw cap gamma sterile	250pcs
LBSST250X	Bottle 250mL PET 25mg sodium thiosulphate spray label gamma sterile	90pcs
LBSST30005X	Container 70mL PP 10mg sodium thiosulphate tablet blue cap label gamma sterile	500pcs
LBSST30132X	Container 120mL PS 10mg sodium thiosulphate labelled blue PE screw cap gamma sterile	300pcs
LBSST32252X	Container 250mL PS 25mg sodium thiosulphate blue cap round gamma sterile	140pcs
LBSST500XS	Bottle 500mL PET 50mg sodium thiosulphate sprayed white screw cap gamma sterile	72pcs
LBS33036THIO	Container 400mL PET square blue wadded cap sodium thiosulphate gamma sterile	60pcs



Coverslips

Item	Description	Unit
MENCS22221GP	Coverslip 22x22mm No.1	1000pcs
MENCS22401GP	Coverslip 22x40mm No.1	1000pcs
MENCS24501GP	Coverslip 24x50mm No.1	1000pcs
MENCS24601GP	Coverslip 24x60mm No.1	1000pcs



Culture Tubes

Kimble have a range of disposable and reusable culture tubes constructed from resilient borosilicate glass. Choose from the reusable culture tubes with PP screw cap and PTFE faced rubber liner, or the rimless disposable culture tubes.

Item	Description	Unit
KIM45066A-13100	Culture tube reusable screw cap teflon liner 13x100mm	72pcs
KIM45066A-16100	Culture tube reusable screw cap teflon liner 16x100mm	72pcs
KIM45066A-16125	Culture tube reusable screw cap teflon liner 16x125mm	72pcs
KIM45066A-16150	Culture tube reusable screw cap teflon liner 16x150mm	72pcs
KIM45066A-20125	Culture tube reusable screw cap teflon liner 20x125mm	48pcs
KIM45066A-25150	Culture tube reusable screw cap teflon liner 25x150mm	36pcs
KIM45066A-25200	Culture tube reusable screw cap teflon liner 25x200mm	36pcs
CSE73500-1075	Tube culture 10x75mm rimless disposable	1000pcs
CSE73500-1275	Tube culture 12x75mm rimless disposable	1000pcs
CSE73500-13100	Tube culture 13x100mm rimless disposable	1000pcs



Consumables & General Labware

Culture Tubes (Continued)

Item	Description	Unit
CSE73500-16100	Tube culture 16x100mm rimless disposable	1000pcs
CSE73500-16125	Tube culture 16x125mm rimless disposable	1000pcs
CSE73500-16150	Tube culture 16x150mm rimless disposable	1000pcs
CSE73500-18150	Tube culture 18x150mm rimless disposable	500pcs
CSE73500-25150	Tube culture 25x150mm rimless disposable	500pcs



Cuvettes

Use these economical, polystyrene, single-use cuvettes for visible spectral analysis.

Item	Description	Unit
LPI111117	2-4mL standard 10mm path length 340-800nm disposable PS	100pcs
LPI112117	0.5-2mL semi-micro 10mm path length 340-800nm disposable PS	100pcs
SIE1939	4.5mL standard PMMA 10mm path length 285-750nm	100pcs
SIE1948	1.5mL semi-micro PMMA 10mm path length 285-750nm	100pcs
STA1/G/10/M/2	Starna Type 1 Glass 10mm path length 3.5mL cuvette matched pair	EA
STA1/G/10	Starna Type 1 Glass 10mm path length 3.5mL cuvette	EA
STA1/Q/2/M/2	Starna Type 1 Quartz 2mm path length 0.700mL cuvette matched pair	EA



Draining Racks

Item	Description	Unit
MET4050	Rack draining 570x600mm wall mounted white nylon coated	EA
MET4055	Rack draining 360x500mm bench mounted white nylon coated	EA
SIE213	Draining rack PS 72 peg capacity 450x630x110mm	EA



Flasks (Erlenmeyer)

Kimble Erlenmeyer flasks are manufactured from borosilicate glass and are graduated. Produced to ISO 1773 standards.

Item	Description	Unit
KIM26500-50	Flask Erlenmeyer 50mL narrow neck graduated	EA
KIM26500-125	Flask Erlenmeyer 125mL narrow neck graduated	EA
KIM26500-250	Flask Erlenmeyer 250mL narrow neck graduated	EA
KIM26500-300	Flask Erlenmeyer 300mL narrow neck graduated	EA
KIM26500-500	Flask Erlenmeyer 500mL narrow neck graduated	EA
KIM26500-1000	Flask Erlenmeyer 1000mL narrow neck graduated	EA
KIM26500-2000	Flask Erlenmeyer 2000mL narrow neck graduated	EA



Consumables & General Labware

Flasks (Volumetric)

Item	Description	Unit
FSBFB50141	Flask volumetric 5mL class A 10/19 PE stopper	EA
FSBFB50143	Flask volumetric 10mL class A 10/19 PE stopper	EA
FSBFB50145	Flask volumetric 20mL class A 10/19 PE stopper	EA
FSBFB50147	Flask volumetric 25mL class A 10/19 PE stopper	EA
FSBFB50149	Flask volumetric 50mL class A 12/21 PE stopper	EA
FSBFB50151	Flask volumetric 100mL class A 12/22 PE stopper	EA
FSBFB50153	Flask volumetric 200mL class A 14/23 PE stopper	EA
FSBFB50155	Flask volumetric 250mL class A 14/23 PE stopper	EA
FSBFB50157	Flask volumetric 500mL class A 19/26 PE stopper	EA
FSBFB50159	Flask volumetric 1000mL class A 24/29 PE stopper	EA
FSBFB50161	Flask volumetric 2000mL class A 29/32 PE stopper	EA



Funnels (Glass)

Item	Description	Unit
PYR2140/02M	Funnel filter 45mm short stem, suits filter paper 55-70mm	EA
PYR2140/03M	Funnel filter 70mm short stem, suits filter paper 70-90mm	EA
PYR2140/10M	Funnel filter 105mm short stem, suits filter paper 150-185mm	EA
PYR2140/12M	Funnel filter 150mm short stem, suits filter paper 240-270mm	EA



Funnels (Separating)

Item	Description	Unit
NAL4300-0500	Funnel separating 500ml PP with TFE stopcock 65mm stem length	4pcs
NAL4300-1000	Funnel separating 1000ml PP with TFE stopcock 65mm stem length	4pcs



Gloves

Item	Description	Unit
TFS148060	Gloves latex powder free x-small 100/box	EA
TFS148070	Gloves latex powder free small 100/box	EA
TFS148080	Gloves latex powder free medium 100/box	EA
TFS148090	Gloves latex powder free large 100/box	EA
TFS148095	Gloves latex powder free x-large 100/box	EA
TFS173060	Gloves nitrile powder free x-small 200/box	EA
TFS173070	Gloves nitrile powder free small 200/box	EA
TFS173080	Gloves nitrile powder free medium 200/box	EA
TFS173090	Gloves nitrile powder free large 200/box	EA
TFS173095	Gloves nitrile powder free x-large 200/box	EA



Consumables & General Labware

Hydrometers

Item	Description	Unit
CARNAL010	Alcoholometer 0-10 x 0.1%	EA
CARNAL020	Alcoholometer 10-20 x 0.2%	EA
CARNAL2040	Alcoholometer 20-40 x 0.2%	EA
CARNAL0100	Alcoholometer 0-100 x 2%	EA
CARBBX010	Brix hydrometer 0-10 x 0.1Bx	EA
CARBBX1020	Brix hydrometer 10-20 x 0.1Bx	EA
CARBBX2030	Brix hydrometer 20-30 x 0.1Bx	EA
CARBBX030	Brix hydrometer 0-30 x 0.5Bx	EA
CAROBM210	Baume Hydrometer -2 to 10 x 0.1Be	EA
CAROBM1020	Baume Hydrometer 10 to 20 x 0.1Be	EA
CAROBM030	Baume Hydrometer 0 to 30 x 0.5Be	EA



Hydrometer Jars

Item	Description	Unit
NAL6230-0500	Hydrometer jar polymethylpentene 50.8 x 390mm	6pcs
CARHJ40	Hydrometer jar glass ungraduated 300mm x 40mm	EA

Inoculating Loops, Needles & Spreaders

Item	Description	Unit
LBSP7201	Inoculating loop PS yellow 1µL sterile 170mm long 20pcs/bag	1000pcs
LBSP7202	Inoculating loop PS blue 10µL sterile 170mm long 20pcs/bag	1000pcs
LBSSPR001X	Spreader PS blue sterile 25pcs/sleeve	3000pcs
LBSSPR002X	Spreader T bar PS blue sterile 10pcs/sleeve	2000pcs

Sterile disposable culture loops - shelf life 3 years

TS/1-A	Inoculating loop 10µL white flexible sterile 20pcs/bag	1000pcs
TS/1-C	Inoculating loop 10µL white flexible sterile individually wrapped	1000pcs
TS/2-A	Inoculating loop 10µL green firm sterile 20pcs/bag	1000pcs
TS/2-C	Inoculating loop 10µL green firm sterile individually wrapped	1000pcs
TS/30-D	Inoculating loop 5µL red flexible sterile 20pcs/bag	1000pcs
TS/30-G	Inoculating loop 5µL red flexible sterile individually wrapped	1000pcs
TS/4-A	Inoculating loop 1µL clear firm sterile 20pcs/bag	1000pcs
TS/4-C	Inoculating loop 1µL clear firm sterile individually wrapped	1000pcs
TS/4-D	Inoculating loop 1µL yellow flexible sterile 20pcs/bag	1000pcs
TS/4-G	Inoculating loop 1µL yellow flexible sterile individually wrapped	1000pcs



Consumables & General Labware

Inoculating Loops, Needles & Spreaders (Continued)

Item	Description	Unit
Sterile disposable inoculating needles - shelf life 2 years		
TS/3-A	Inoculating needle clear firm sterile 20pcs/bag	1000pcs
TS/3-C	Inoculating needle clear firm sterile individually wrapped	1000pcs
Sterile disposable spreaders - shelf life 3 years		
TS/30-A	Spreader L shaped clear firm sterile 5pcs/bag	1000pcs
TS/30-B	Spreader L shaped clear firm sterile individually wrapped	1000pcs
TS/30-C	Spreader L shaped clear firm sterile 25pcs/bag	1000pcs
TS/30-ES	Spreader L shaped blue flexible sterile 25pcs/bag	1000pcs



Jugs

Strong and stable, these high clarity polypropylene autoclavable jugs have well-formed pouring spouts and handles with clearly defined graduations.

Item	Description	Unit
LBSJUG250	Jug 250mL PP with printed graduations autoclavable	EA
LBSJUG500	Jug 500mL PP with printed graduations autoclavable	EA
LBSJUG1000	Jug 1L PP with printed graduations autoclavable	EA
LBSJUG2000	Jug 2L PP with printed graduations autoclavable	EA
LBSJUG3000	Jug 3L PP with printed graduations autoclavable	EA
LBSJUG5000	Jug 5L PP with printed graduations autoclavable	EA



Label Printer

Item	Description	Unit
BAD139538	BMP21 LAB Printer	EA
BAD110932	Self Laminating Polyester Label, Black on White, 12.7mm x 6.4m	EA
BAD110935	High Performance Polyester Label, Black on White, 12.7mm x 6.4m	EA
BAD110894	Nylon Cloth Label, Black on White, 12.7mm x 6.4m	EA
BAD110900	Clear Polyester Label, Black on Clear, 12.7mm x 6.4m	EA
BAD139747	Clear Polyester Label, White on Clear, 12.7mm x 6.4m	EA



Magnetic Stirrer Bars

Item	Description	Unit
COW1.513.3	PTFE Stirrer Bar, Octahedral 13mm x 3mm	EA
COW001.513	PTFE Stirrer Bar, Octahedral 13mm x 8mm	EA
COW001.513.10	PTFE Stirrer Bar, Octahedral 13mm x 10mm	EA
COW1.515	PTFE Stirrer Bar, Octahedral 15mm x 8mm	EA
COW1.525	PTFE Stirrer Bar, Octahedral 25mm x 8mm	EA
COW001.525.10	PTFE Stirrer Bar, Octahedral 25mm x 10mm	EA
COW001.528	PTFE Stirrer Bar, Octahedral 28mm x 8mm	EA
COW001.538	PTFE Stirrer Bar, Octahedral 38mm x 8mm	EA
BAT37110.0178	Spinbar Teflon PTFE Octagon 47.6mm x 9.5mm	EA

Consumables & General Labware

Item	Description	Unit
COW001.538.10	PTFE Stirrer Bar, Octahedral 38mm x 10mm	EA
COW001.551.10	PTFE Stirrer Bar, Octahedral 51mm x 10mm	EA
COW001.564	PTFE Stirrer Bar, Octahedral 64mm x 10mm	EA
COW001.575	PTFE Stirrer Bar, Octahedral 75mm x 13mm	EA



Markham Apparatus

Item	Description	Unit
CARMS3	Markham still body with stopper, S29 ball joint	EA
CARSHC	Markham splash head cyclone with 24/29 cone, S29 cup	EA
CARC7/5/200	Condenser, inland revenue with 24/29 joints, 200mm effective length	EA
CARAD11/5	Markham adapter, delivery straight with 24/29 joint	EA
CARCLK24	Clip joint, green plastic, suits 24/29 joints (2 required)	EA
CARCLS29	Clip, spherical joint, S29	EA
CARMA-1	Markham apparatus set complete includes still body, splash-head, condensor, adaptors and clips	EA



Rankine Apparatus

Item	Description	Unit
CARF1/5/100	Rankine round bottom flask with 24/29 neck, 100mL	EA
CARAD6/55.14	Rankine adapter swan neck with 24/29 joints, GL14 screwthread	EA
CART2/250	Rankine air inlet tube 7mm x 250mm tapered end	EA
CARC7/5/200	Condenser, inland revenue with 24/29 joints, 200mm effective length	EA
CARAD15.5.13	Rankine adapter bridge S13 ball to 24/29 cone	EA
CARCLS13	Clip, spherical joint, S13 steel with adjuster	EA
CARBUB1/3.13	Rankine Bubbler with 14/23 cone, S13 cup	EA
CARF7/3/50	Rankine flask pear shape 2 necks with 14/23 sockets, 50mL	EA
CARCLK24	Clip joint, green plastic, suits 24/29 joints (3 required)	EA
CARCLK14	Clip joint, yellow plastic, suits 14/23 joints (2 required)	EA
CARSETRA1	Rankine apparatus set complete includes flasks, adaptors, condensers, bubblers and all clips	EA



Measuring Cylinders (Glass)

Class A with white graduations.

Item	Description	Unit
PYR3201/04	Cylinder measuring 10mL	EA
PYR3201/06	Cylinder measuring 25mL	EA
PYR3201/08	Cylinder measuring 50mL	EA
PYR3201/10	Cylinder measuring 100mL	EA
PYR3201/12	Cylinder measuring 250mL	EA
PYR3201/14	Cylinder measuring 500mL	EA

Consumables & General Labware

Item	Description	Unit
PYR3201/16	Cylinder measuring 1000mL	EA
PYR3201/18	Cylinder measuring 2000mL	EA

Measuring Cylinders (Plastic)

With moulded graduations hexagonal base.

Item	Description	Unit
LBS80303H-03	Cylinder measuring 50mL PP	EA
LBS80303H-04	Cylinder measuring 100mL PP	EA
LBS80303H-05	Cylinder measuring 250mL PP	EA
LBS80303H-06	Cylinder measuring 500mL PP	EA
LBS80303H-07	Cylinder measuring 1000mL PP	EA
NAL3664-0050	Cylinder measuring 50ml PP, graduated	EA
NAL3664-0100	Cylinder measuring 100ml PP, graduated	EA
NAL3664-0250	Cylinder measuring 250ml PP, graduated	EA
NAL3664-0500	Cylinder measuring 500ml PP, graduated	EA
NAL3664-1000	Cylinder measuring 1000ml PP, graduated	EA



Parafilm

Item	Description	Unit
PAAPM996	M Sealing Film 4inch x 125ft roll	EA
PAAPM992	M Sealing Film 2inch x 250ft roll	EA

Petri Dishes

Our LabServ manufacturing plant produces a selection of good quality, disposable, polystyrene petri dishes. Other size sterility and configuration options also available.

Item	Description	Unit
LBS60001	Dish petri PS 90x15mm manual vented with stacking rim sleeves of 20	500pcs
LBS60001X	Dish petri PS 90x15mm manual sterile vented with stacking rim sleeves of 20	500pcs



Pipettes (Bulb/Volumetric)

One mark class AS amber stain.

Item	Description	Unit
PYRPRA720AS	Pipette volumetric 1mL	2pcs
PYRPRA722AS	Pipette volumetric 2mL	2pcs
PYRPRA724AS	Pipette volumetric 5mL	2pcs
PYRPRA726AS	Pipette volumetric 10mL	2pcs



Consumables & General Labware

Item	Description	Unit
PYRPRA728AS	Pipette volumetric 15mL	2pcs
PYRPRA732AS	Pipette volumetric 25mL	2pcs
PYRPRA734AS	Pipette volumetric 50mL	2pcs
PYRPRA736AS	Pipette volumetric 100mL	2pcs

Pipettes (Graduated)

Type 2 class AS amber graduations.

Item	Description	Unit
PYRPRA810AS	Pipette graduated 1mL	2pcs
PYRPRA812AS	Pipette graduated 2mL	2pcs
PYRPRA814AS	Pipette graduated 5mL	2pcs
PYRPRA816AS	Pipette graduated 10mL	2pcs
PYRPRA818AS	Pipette graduated 25mL	2pcs

Pipettes (Serological)

LabServ offers a wide range of bulk serological pipettes for liquid transfer. LabServ disposable serological pipettes are of sterile optically clear polystyrene construction. They are colour coded with bidirectional graduations for flexibility of the user. Wrapped with lot number traceability.

Item	Description	Unit
LBSCN857101	1mL PS individually wrapped sterile	800pcs
LBSCN857102	2mL PS individually wrapped sterile	600pcs
LBSCN857105	5mL PS individually wrapped sterile	400pcs
LBSCN857110	10mL PS individually wrapped sterile	400pcs
LBSCN857125	25mL PS individually wrapped sterile	200pcs
LBSCN857201	1mL PS sterile bulk packs of 25pcs per pack	1000pcs
LBSCN857202	2mL PS sterile bulk packs of 25pcs per pack	800pcs
LBSCN857205	5mL PS sterile bulk packs of 25pcs per pack	600pcs
LBSCN857210	10mL PS sterile bulk packs of 25pcs per pack	600pcs
LBSCN857225	25mL PS sterile bulk packs of 25pcs per pack	400pcs
NUN170358N	50ml PS individually wrapped paper pack, PET filter, sterile	100pcs
NUN170369N	50ml PS individually wrapped plastic pack, PET filter, sterile	100pcs
NUN170376N	50ml PS bulk pack of 25pcs per pack, PET filter, sterile	100pcs

Pipette Controllers

Includes table stand, wall mount, 1mL pipette stabiliser, universal charging transformer.

Item	Description	Unit
MAT9511	S1 cordless serological pipettor clear	EA
MAT9521	S1 cordless serological pipettor blue	EA



Consumables & General Labware

Item	Description	Unit
MAT9531	S1 cordless serological pipettor red	EA
MAT9541	S1 cordless serological pipettor green	EA

Pipettes (Transfer)

Prevent cross-contamination with Thermo Scientific™ Samco™ General-Purpose Transfer Pipettes, which have a built-in bulb that eliminates the need for plugging, packaging and autoclaving. Samco transfer pipettes are made of shatterproof, non-toxic polyethylene. Their low-affinity surface reduces loss of cells and valuable proteins due to binding. They can be frozen in liquid nitrogen, or heat sealed for storage.

Item	Description	Unit
BER222	Transfer pipette graduated 1mL large bulb	500pcs
BER222-1S	Transfer pipette graduated 1mL sterile individually wrapped	500pcs
BER225	Transfer pipette graduated 3mL large bulb	500pcs
BER225-1S	Transfer pipette graduated 3mL sterile individually wrapped	500pcs

Pipettes

We offer a complete range of Finnpiquette® manual and electronic pipettors. Thermo Scientific Finnpiquette F1 combine all the features that Finnpiquettes are famous for plus a patented soft-touch tip ejection, an adjustable finger rest and an antimicrobial surface to prevent contamination of samples and pipettors. Single and multichannel fixed volume and variable volume options available.

Item	Description	Unit
FNP4642060	Finnpiquette F2 2-20µL manual pipettor	EA
FNP4642080	Finnpiquette F2 20-200µL manual pipettor	EA
FNP4642090	Finnpiquette F2 100-1000µL manual pipettor	EA
FNP4642100	Finnpiquette F2 0.5-5mL manual pipettor	EA
FNP4642110	Finnpiquette F2 1-10mL manual pipettor	EA
FNP4641180	Finnpiquette F1 ClipTip 2-20µL manual pipettor	EA
FNP4641210	Finnpiquette F1 ClipTip 20-200µL manual pipettor	EA
FNP4641230	Finnpiquette F1 ClipTip 100-1000µL manual pipettor	EA
FNP4661110	Finnpiquette F1 ClipTip 8-channel 1-10µL manual pipettor	EA
FNP4661120	Finnpiquette F1 ClipTip 8-channel 5-50µL manual pipettor	EA
FNP4661130	Finnpiquette F1 ClipTip 8-channel 10-100µL manual pipettor	EA
FNP4661140	Finnpiquette F1 ClipTip 8-channel 30-300µL manual pipettor	EA

Pipette Tips

Item	Description	Unit
FNP94300120	Pipette tip yellow 0.2-200µL bulk	1000pcs
FNP94300220	Pipette tip blue 100-1000µL bulk	1000pcs
FNP9400260	Finnpip 250 0.5-250µL natural racked, 10 racks of 96 tips	960pcs



Consumables & General Labware

Item	Description	Unit
FNP9401110	Finntip 1000 100-1000µL natural racked, 10 racks of 96 tips	960pcs
QSP095	Pipette tip 1-5mL natural bulk	250pcs
FNP9402151	Finntip 10mL 1000-10000µL natural	100pcs
FNP94410210	ClipTip 20 non sterile, 10 racks of 96 tips	960pcs
FNP94410217	ClipTip 20 Reload non sterile, 10 inserts of 96 tips	960pcs
FNP94410250	ClipTip 50 non sterile, 10 racks of 96 tips	960pcs
FNP94410257	ClipTip 50 Reload non sterile, 10 inserts of 96 tips	960pcs
FNP94410310	ClipTip 200 non sterile, 10 racks of 96 tips	960pcs
FNP94410317	ClipTip 200 Reload non sterile, 10 inserts of 96 tips	960pcs
FNP94410510	ClipTip 300 non sterile, 10 racks of 96 tips	960pcs
FNP94410517	ClipTip 300 Reload non sterile, 10 inserts of 96 tips	960pcs
FNP94410710	ClipTip 1000 non sterile, 10 racks of 96 tips	960pcs
FNP94410717	ClipTip 1000 Reload non sterile, 10 inserts of 96 tips	960pcs



Safety Glasses

Item	Description	Unit
BOL1652301	Bollé Rush Safety Glasses Clear Lens	EA
BOL1667201	Bollé Bandido Safety Glasses Clear Lens	EA
BOL1668201	Bollé Blade Safety Glasses Clear Lens	EA
BOL1650515	Bollé Override Cover Spectacles Clear Lens	EA
BOL1661100	Bollé Extreme Safety Glasses With Sliver Frames & Clear Lenses	CS10



Slides

Item	Description	Unit
FNN7101WT	Slide plain ground white edge 1-1.2mm 26x76mm	2500pcs
FNN7105WT	Slide single frost ground edge white 26x76mm	2500pcs

Spatulas

Item	Description	Unit
MET2715-160	Micro spoon with spatula 150mm long	EA
MET2712-150	Micro spoon with spatula 150mm long	EA
N-E9370	Spoon spatula 160mm	10pcs
FSH14-648-231	Digital Spatula Balance scoops and weighs 0 to 300g in one unit	EA



Syringes and Needles

Luer lock options are also available. Contact us for more information.

Item	Description	Unit
TER1NN.1838R	Needle 18g x 38mm	100pcs

Consumables & General Labware

Item	Description	Unit
TER1NN.2138R	Needle 21g x 38mm	100pcs
TER1NN.2325RL	Needle 23g x 25mm	100pcs
TER1NN.2519RA	Needle 25g x 19mm	100pcs
TER1NN.2713R	Needle 27g x 13mm	100pcs
TER6SS.01T	Syringe 1mL leur slip	100pcs
TER6SS.03S	Syringe 3mL leur slip	100pcs
TER6SS.05S	Syringe 5mL leur slip	100pcs
TER6SS.10S	Syringe 10mL leur slip	100pcs
TER6SS.10ES	Syringe 10mL eccentric leur slip	1200pcs
TER6SS.20ES	Syringe 20mL eccentric leur slip	50pcs
TER1SS.30ESE1	Syringe 30mL eccentric leur slip	50pcs
TER1SS.50ESE	Syringe 50mL eccentric leur slip	20pcs



Test Tube Racks

To select the different colours replace the 'X' in code below with
 0 = white, 1 = orange, 2 = yellow, 3 = blue, 4 = green or 5 = red.

Item	Description	Unit
NAL5970-0X13	13mm diameter tube - 72 tubes 6 x 12 array	EA
NAL5970-0X16	16mm diameter tube - 72 tubes 6 x 12 array	EA
NAL5970-0X20	20mm diameter tube - 40 tubes 4 x 10 array	EA
NAL5970-0X25	25mm diameter tube - 40 tubes 4 x 10 array	EA
NAL5970-0X30	30mm diameter tube - 24 tubes 3x 8 array	EA



Timers

Item	Description	Unit
FSH02-261-840	Multi-colored timer four channel count down 100 min with alarm easel bench stand clip and magnet	EA
FSH06-662-45	Lab-top timer count up/down 24 hours with clock for bench use	EA
FSH14-649-14	Timer count up, 60 min count down, 100 min rear mounted magnet pocket clip bench stand LCD display and alarm	EA
FSH14-649-17	Timer four channel count up/down 99 hours clock alarm jumbo LCD display rear mounted magnet clip bench stand	EA

Thermometers (Digital)

Item	Description	Unit
FSH14-648-44	Waterproof lollipop thermometer -50 to 300°C (NIST Traceable)	EA



Consumables & General Labware

Item	Description	Unit
FSH15-077-9E	Waterproof thermometer -50 to 300°C 3m cable with 120mm SS probe (NIST Traceable)	EA
FSH06-664-11	Fridge freezer thermometer -50 to 70°C enclosed with temperature-buffered sensors (NIST Traceable)	EA



Thermometers (General)

76mm immersion 305mm long.

Item	Description	Unit
LBS72/809/0	-10 to 50°C red spirit white back 0.5 divisions	EA
LBS72/808/0	-10 to 110°C red spirit white back 1.0 divisions	EA
LBS72/807/0	-10 to 150°C red spirit white back 1.0 divisions	EA
LBS72/810/0	-10 to 110°C Hg white back 1.0 divisions	EA



Vials and Closures

Item	Description	Unit
THC11090500	Vial Screw 2mL clear no markings	100pcs
THC09150838	Cap Screw ultraclean silicone/PTFE 9mm blue	100pcs
THC18091307	Vial Headspace 20mL screw top clear	100pcs
THC18032063	Cap Screw 18mm Magnetic Sil/PTFE slit septa	100pcs
THC18031414	Cap Screw 18mm Magnetic silicone/PTFE septa	100pcs



Thermo Fisher Scientific is proud to represent Pyrex, Quickfit, MBL & Azlon® for the first time. Trusted by generations of scientists, the original suppliers of the highest quality borosilicate glass and reusable plastics are now available from the World Leader in Serving Science. With over 3000 newly introduced general purpose and specialised lines, Thermo Fisher Scientific is the location to source all your highest quality glassware & plasticware needs with an unsurpassed range from the most recognised names in science.

Download your new 320-page laboratory glassware, plasticware & accessories catalogue now!



thermofisher.com.au/scilabwarecatalogue

Consumables & General Labware

Wash Bottles (Safety)

Thermo Scientific Nalgene labelled LDPE "Right to Know" wash bottles feature the self-venting system that keeps hazardous volatile liquids in the bottle and off the bench.

Item	Description	Unit
NAL2425-0501	Acetone (Red) 500mL	6pcs
NAL2425-0502	Ethyl Alcohol (White) 500mL	6pcs
NAL2425-0503	Methanol (Green) 500mL	6pcs
NAL2425-0504	Isopropanol (Yellow) 500mL	6pcs
NAL2425-0505	Distilled Water (Natural) 500mL	6pcs
NAL2425-0506	Sodium Hypochlorite (White) 500mL	6pcs
AZLWGF336P	Bottle wash 500mL slope shoulder narrow neck LDPE	EA
AZLWGF240P	Bottle wash 1L slope shoulder narrow neck LDPE	EA
AZLBGF130P	Bottle wash 250mL integral narrow neck LDPE	EA
AZLBGF136P	Bottle wash 500mL integral narrow neck LDPE	EA



Wash Bottles - GHS labelled

Thermo Scientific Nalgene "Right to Understand" wash bottles are labelled with new Globally Harmonized System (GHS) icons and cautions.

Item	Description	Unit
NAL2428-0502	500ml Ethanol, LDPE, white closure	6pcs
NAL2428-0503	500ml Methanol, LDPE, green closure	6pcs
NAL2428-0504	500ml Isopropanol, LDPE, yellow closure	6pcs
NAL2428-0505	500ml Distilled Water, LDPE, natural closure	6pcs
NAL2428-0506	500ml Sodium Hypochlorite, LDPE, White closure	6pcs
NAL2428-0507	500ml Deionised Water, LDPE, natural closure	6pcs



Unitary Safety Wash Bottles

Vented unitary bottle design prevents contamination when filling the bottle; no draw tube assembly to pick-up contamination off the bench.

Item	Description	Unit
NAL2436-0501	Acetone (Red) 500mL	EA
NAL2436-0502	Ethyl Alcohol (White) 500mL	EA
NAL2436-0503	Methanol (Green) 500mL	EA
NAL2436-0504	Isopropanol (Yellow) 500mL	EA
NAL2436-0505	Distilled Water (Natural) 500mL	EA
NAL2436-0506	Sodium Hypochlorite (White) 500mL	EA
NAL2436-0507	De-ionised Water (Natural) 500mL	EA



Weigh Boats

Item	Description	Unit
SLNWB30205	Weigh boat standard square micro 44 x 44mm 7mL	500pcs
SLNWB30254	Weigh boat standard square regular 80 x 80mm 100mL	500pcs
SLNWB30304	Weigh boat standard square regular 140 x 140mm 250mL	500pcs



Equipment

Autoclaves	29
Balances	29
Centrifuges	30
CO ₂ monitoring	31
Electrochemistry	32
Homogenisers	33
Hotplates & Hotplate Magnetic Stirrers	34
Incubators	34
Microplate Readers	35
Ovens	35
Oxygen Demand (Biological & Chemical)	36
Platform Rockers Shakers & Mixers	36
Refractometers	37
Refrigerators and Freezers - General Purpose Laboratory	38
Spectrophotometers	38
Turbidimeters	39
Vortex	39
Water Baths	40
Water Purification Systems	41

Absolutely the best quality, most reliable laboratory benchtop equipment from Thermo Scientific and our supply partners.

Autoclaves

The Tuttnauer D-line Benchtop Series of autoclaves are fully automated and microprocessor controlled providing safe, economical and effective sterilisation. A range of sizes available suitable for a variety of applications.

Item Code	Description	Unit
TTN2840ELD	28L Benchtop Autoclave with Built-in Reservoir	EA
TTN3840ELD	52L Benchtop Autoclave with Built-in Reservoir	EA
TTN3870ELD	85L Benchtop Autoclave with Built-in Reservoir	EA



Balances (Sartorius Entris Series)

Get the reliability and quality of a Sartorius balance, but only pay for the features you really need. The new Entris® has been specifically designed to deliver effective and reliable weighing results in your daily work. Offering internal calibration as standard. Other models available on request.


Item	Description	Unit
SARENTRIS224I-1S	Entris Analytical Balance 220g 4dp	EA
SARENTRIS124I-1S	Entris Analytical Balance 120g 4dp	EA
SARENTRIS623I-1S	Entris Precision Balance 620g 3dp	EA
SARENTRIS822I-1S	Entris Precision Balance 820g 2dp	EA
SARENTRIS2201I-1S	Entris Precision Balance 2200g 1dp	EA



Balances (Sartorius Quintix Series)

Sets a new benchmark in every aspect for standard lab balances (19 models available). All units feature internal calibration. Other models available on request.

Item	Description	Unit
Analytical Balances		
SARQUINTIX124-1S	Quintix Analytical Balance 120g 4dp	EA
SARQUINTIX224-1S	Quintix Analytical Balance 220g 4dp	EA
Precision Balances		
SARQUINTIX513-1S	Quintix Precision Balance 510g 3dp	EA
SARQUINTIX612-1S	Quintix Precision Balance 610g 2dp	EA



Centrifuges

Thermo Fisher Scientific centrifuge systems deliver outstanding performance and reliability in the lab. From microcentrifuges and benchtop instruments, to advanced floor models and innovative carbon fiber rotors - all offer exceptional quality for consistent results.

Item	Description	Unit
THR75002410	Heraeus Pico 17 microcentrifuge, includes 24 x 1.5/2.0mL rotor with ClickSeal biocontainment lid, max speed 17000 x g	EA
THR75007211	Heraeus Megafuge 8 centrifuge, requires rotor and adapters	EA
THR75004210	Heraeus Multifuge X1 centrifuge, requires rotor and adapters	EA
THR75004500	Heraeus Multifuge X3 centrifuge, requires rotor and adapters	EA



Microcentrifuges

Accelerate your micro-volume separations with our compact, safe, and easy-to-use microcentrifuges. From intuitive controls, to innovative features our microcentrifuges deliver the perfect blend of capabilities to support your micro-volume preparation – all in a remarkably small footprint.

Item	Description	Unit
THR75004061	mySpin 6 ventilated microcentrifuge with 6x 1.5/2.0mL, 16x 0.2mL rotor & adapters Max RCF: 2,000 x g	EA
THR75004081	mySpin 12 ventilated microcentrifuge with 12x 1.5/2.0mL, 32x 0.2mL rotor & adapters Max RCF: 9,800 x g	EA
THR75002410	Heraeus Pico 17 ventilated microcentrifuge with 24 x 1.5/2.0mL rotor & biocontainment lid Max RCF: 17000 x g	EA
THR75002420	Heraeus Fresco 17 refrigerated microcentrifuge with 24 x 1.5/2.0mL rotor & biocontainment lid Max RCF: 17000 x g	EA



Small Bench Centrifuges

The Thermo Scientific™ Heraeus Megafuge 8 series centrifuges provide exceptional capacity in a compact design with a smart, simple interface and the flexibility to support a range of applications. Featuring secure, push-button Auto-Lock™ rotor exchange, one-touch operation with pre-saved protocols and a unique 8x 15/50ml swinging bucket rotor, the Megafuge 8 series is the perfect fit for your laboratory.

Item	Description	Unit
THR75007211	Megafuge 8 Ventilated Small Bench Centrifuge	EA
THR75007214	Megafuge 8R Refrigerated Small Bench Centrifuge	EA
THR75005761	High Capacity Cell Culture Rotor Package 8x 15mL or 8x 50mL Conical Tubes	EA



Safety Facts

Why is monitoring CO₂ critical?

Exposure to high levels of Carbon Dioxide (CO₂) can be dangerous to workers in Wineries. Care must be taken to protect workers from over exposure. CO₂ is twice as heavy as air; therefore it will sink to the bottom of a room forming potentially deadly pools of gas that will displace oxygen. Areas of concern within Wineries include pits, sumps and storage tanks, as well as fermentation rooms, barrel cellars and bottling rooms.

How the Sensepoint XCD can help

The detection of an excessive level of CO₂ in the monitored area produces an alarm, alerting staff working in the affected area and initiating a response from mechanical ventilation systems.

Advantages

- Local back-lit LCD display of the CO₂ concentration in the area.
- Provides tri-colored alarm indication – red backlight alerts nearby staff to an immediate hazard
- A variety of output signal protocols (Analogue, MODBUS, Alarm Relays.)



Item	Description	Unit
SGRSPXCDALMB1	MD200 COD photometer	EA

To find out more, contact 1300 735 295

Electrochemistry (Conductivity/TDS)

Choose from a wide range of conductivity meters and kits. Handheld and bench top testers available. Contact us for electrode buffers and storage solutions.

Item	Description	Unit
Testers		
EUTELITECTSPIN	Thermo Scientific Eutech Elite conductivity/TDS/salinity pocket tester with pin sensor	EA
Portable		
EUT01X651633	CON 450 Conductivity/TDS/Salinity waterproof meter with K=1.0 electrode, solutions, and carrying case.	EA
ORISTARA2225W	Orion Star A222 portable conductivity meter kit with 013010MD conductivity DuraProbe solutions and carry case.	EA
Bench Top		
EUT01X543905W	CON 2700 Conductivity/TDS/Temperature meter with conductivity probe and integral probe holder.	EA
ORISTARA2125W	Orion Star A212 benchtop conductivity meter kit with 013005MD conductivity DuraProbe electrode arm and 1413µS/cm standard.	EA



Electrochemistry (pH)

If you need more than a pH test strip, then choose from the handheld and benchtop pH testers listed below. Contact us for electrodes buffers and storage solutions.

Item	Description	Unit
pH Testers		
EUTELITEPH	Thermo Scientific Eutech Elite pH pocket tester	EA
Handheld		
EUT01X651630	Eutech pH 450 portable pH meter kit with double junction pH electrode, temperature probe, buffers and carry case.	EA
ORISTARA2215W	Orion Star A221 portable pH meter kit with 8107UWMMD ROSS pH/ATC triode protective armour holder hard field case and solutions	EA
Bench Top		
EUT01X53902W	Eutech pH 2700 benchtop pH meter with double junction pH probe temperature probe and electrode holder.	EA
ORISTARA2115W	Orion Star A211 benchtop pH meter kit with 8320BNUMD ROSS pH/ATC triode buffers, storage solution, cleaning solution and electrode stand.	EA



Electrochemistry – (Thermo Scientific Orion AquaMate)

Item	Description	Unit
ORIAQ8000	Orion AquaMate 8000 UV/visible spectrophotometer with preprogrammed methods for water analysis on USB	EA

Electrochemistry – (Portable Dissolved Oxygen Meter)

These portable Dissolved Oxygen meters offer the advantages of an optical probe - no membrane, no electrolyte and are not affected by sulphides, resulting in less calibration and maintenance.

Item	Description	Unit
ORISTARA3235W	Orion Star A323 portable Dissolved Oxygen meter with RDO probe with 3m cable and probe guard. Includes calibration sleeve, batteries, and carry case.	EA
AQRDO-200/ AQUAPLUS	AquaRead DO-200 meter with AquaPlus Optical Dissolved Oxygen Probe. System also includes barometric pressure sensor. Integral conductivity and temperature. Data logging. Standard 3m cable. Extension cables are available. Removable metal shield fits around the probe. Minus 5 to 50°C probe operating temperature.	EA



Electrochemistry – (Multiparameter)

Test multiple parameters with these portable and benchtop meters below.

Item	Description	Unit
Tester		
EUTELITEPCTS	Thermo Scientific Eutech Elite pH/conductivity/TDS/salinity pocket tester	EA
Portable		
EUT01X430902W	Eutech CyberScan PCD650 Portable pH/Conductivity/Dissolved Oxygen Meter. All parameters can be displayed on LCD screen. Wide measurement ranges. Logs up to 500 readings.	EA
ORISTARA3295W	Orion Star A329 portable pH/ISE/cond/RDO meter kit with pH RDO and conductivity probes armour hard field case and solutions.	EA
Bench Top		
ORIVSTAR92W	Orion VERSA STAR benchtop pH/ISE/cond/DO kit with pH, DO and conductivity probes and electrode stand.	EA



Homogenisers

Handheld and benchtop models available with a selection of probes and stands to fit your application.

Item Code	Description	Unit
BRTEQ02520-300-RD000.0	Precellys Evolution Tissue Homogeniser with 24x2mL/0.5mL holder	EA
OMNTHP220	Homogeniser package includes 3 probes & stand	EA



Hotplates & Hotplate Magnetic Stirrers

Available in three surface sizes the Cimarec+ series of digital stirrers, hot plates and hot plate stirrers offer precise stirring control, exceptional safety and superior temperature for your routine protocols.

Features:

- Ceramic top for excellent chemical resistance
- Unique HOT TOP warning system protects from accidental burning
- Thermo Scientific Stir Trac technology for low speed and strong magnetic coupling
- Maximum temperature 540°C



Item	Description	Unit
BARSP88854105	Cimarec+ stirrer hotplate small ceramic top 10.8 x 10.8 cm	EA
BARSP88857106	Cimarec+ stirrer hotplate medium ceramic top 18.4 x 18.4 cm	EA

Incubators (Benchtop)

Heratherm Microbiological Incubators are available in 60, 100 or 180 litre with 2 airflow technologies (gravity or our new dual convection). They have an inner glass door, small foot print, temperature range ambient +5°C up to 105°C. Also available compact 18 litre model for work-around ambient temperatures. Other models available on request.



Item	Description	Unit
HEA50125882	IMC18 HERATHERM Compact Incubator mechanical convection 17-40°C 18L	EA
HEA51028130	IGS60 HERATHERM General Protocol Incubator Gravity Convection 75L ambient +5°C to 75°C	EA
HEA51028131	IGS100 HERATHERM General Protocol Incubator Gravity Convection 117L ambient +5°C to 75°C	EA
HEA51028132	IGS180 HERATHERM General Protocol Incubator Gravity Convection 194L ambient +5°C to 75°C	EA

Incubators (Large Capacity)

Heratherm Large Capacity Microbiological Incubators are available in 400 or 750 litre size with 2 airflow technologies (gravity or mechanical convection). Inner glass door. Temperature range ambient+5°C up to 105°C.

Item	Description	Unit
HEA51029322	IGS400 HERATHERM General Protocol Incubator Gravity Convection 400L ambient +5°C to 75°C	EA
HEA51029334	IGS750 HERATHERM General Protocol Incubator Gravity Convection 750L ambient +5°C to 75°C	EA
HEA51029325	IMH400-S HERATHERM Advanced Security Protocol Incubator Mechanical Convection 400L ambient +5°C to 105°C	EA
HEA51029337	IMH750-S HERATHERM Advanced Security Protocol Incubator Mechanical Convection 750L ambient +5°C to 105°C	EA



Incubators (Refrigerated)

Heratherm Refrigerated Incubators have a temperature range of +5 to +70°C. With peltier technology for energy savings and reliable stable temperature control.

Item Code	Description	Unit
HEA51031562	IMP180 Heratherm refrigerated incubator 178L Temp range +5 to 70°C	EA
HEA51031565	IMP400 Heratherm refrigerated incubator 381L Temp range +5 to 70°C	EA



Microplate Readers (Photometric and Spectrophotometric)

Thermo Fisher Scientific offers a full range of microplate readers suitable for high throughput sample analysis of wine and grape juice as well as stability tests including glucose, malic acid, L-lactic acid and acetic acid. Microplate readers are capable of testing 96 simultaneous samples in a microplate format. Both filter detection and scanning detection models available.

Thermo Scientific Multiskan GO

The Thermo Scientific Multiskan GO UV/Vis microplate spectrophotometer is ideal for labs doing multiple assays. It offers free wavelength selection (200nm to 1000nm) and reads microplates and cuvettes.



Item Code	Description	Unit
MIB51119200	Multiskan GO without Cuvette	EA
MIB51119300	Multiskan GO with Cuvette	EA

Thermo Scientific Multiskan FC

The Multiskan FC microplate photometer is a reliable and robust instrument for a wide variety of applications. It is a filter based system capable of detection in the wavelength range of 340 to 850nm and reads 96 well plates in a quick 5 seconds.

Item Code	Description	Unit
MIB51119000	Multiskan FC	EA
MIB51119100	Multiskan FC with incubator	EA

Ovens (Benchtop)

Heratherm Ovens are available in three different models: Each unit is available in 60, 100 and 180 litre sizes with a choice of gravity or mechanical convection. Temperature range of 50°C up to 330°C. Other models available on request.

Item	Description	Unit
HEA51028148	OMS60 HERATHERM General Protocol Mechanical Convection 66L 50°C to 250°C	EA
HEA51028149	OMS100 HERATHERM General Protocol Mechanical Convection 104L 50°C to 250°C	EA
HEA51028150	OMS180 HERATHERM General Protocol Mechanical Convection 179L 50°C to 250°C	EA

Oxygen Demand, Chemical

The chemical oxygen demand test measures the equivalent oxygen to the amount of organic matter oxidisable by potassium dichromate in a 50% sulphuric acid solution. It is used as a quick means to estimate BOD for process control.

Item	Description	Unit
TIN2892502	MD200 COD photometer	EA
TIN2418940	RD 125 COD reactor with 24x15mm holes and integrated protection cover	EA
TSPCODL00B150	Low range COD vials 3- 150 mg/L pk/150	EA
TSPCODH00B150	Medium range COD vials 20 - 1500 mg/L pk/150	EA
TSPCODHP0B150	High range COD cials 200 - 15000 mg/L pk/150	EA

Oxygen Demand, Biological

Biological oxygen demand (BOD) is a measure of the amount of oxygen consumed by a wastewater discharge. Required for reporting purposes for waste discharges.

Item	Description	Unit
TIN2444460	BD600 respirometric BOD system includes flasks, stirring bars, overflow bottles, gaskets, and nitrification inhibitor. Features remote control. Requires mains power and to be placed inside and incubator.	EA

Platform Rockers Shakers & Mixers

Whether you require simple rocking or orbital mixers Thermo Fisher Scientific offer a wide choice of devices to suit your needs.

Item	Description	Unit
BAR88880026	Orbital shaker 50-300rpm compact digital, requires platform and clamps	EA
BAR88880020	Compact digital rocker 5-100rpm with timer tilt angle adjustable 7°-13° tray 27.6 x 27.6cm	EA
BAR88880022	Compact digital wave rotator 8-100rpm with timer tilt angle adjustable 0°-13° tray 27.6 x 27.6cm	EA



Refractometers (Digital Handheld)

Optical refractometers continue to be a quick and easy way to check Brix, but why not upgrade to a handheld digital unit with our range from market leader Bellingham+Stanley.

Each unit offers temperature compensation, quick and simple in-situ calibration as well as a robust design with replaceable batteries for ease of use. Small and light, they are ideal for factory or field-based testing.

Item	Description	Unit
BLS38-02	Direct reading digital hand held refractometer. Temperature compensated. Display Brix 0-54%	EA
BLS38-A1	OPTi Brix 95+ refractometer; 0 - 95% Brix	EA
BLS38-41	OPTi Wine AP Refractometer % Mass w/w (0-35) and Alcohol Probable (AP) (0-22) ranges	EA
BLS38-42	OPTi Wine OE-D Refractometer % Mass w/w (0-35) and Oechsle (German) (30-130) ranges	EA
BLS38-43	OPTi Wine OE-CH Refractometer % Mass w/w 0-35 and Oechsle (Swiss) (0-130) ranges	EA
BLS38-44	OPTi Wine KMW Refractometer % Mass w/w 0-35 and KMW (Babo) (0-25) ranges	EA
BLS38-45	OPTi Duo Wine Baume refractometer % Mass w/w 0-35 and °Baumé (0-28) ranges	EA
BLS38-46	OPTi Wine Baume/AP % Mass w/w 0-28 and Alcohol Probable (AP) (0-22) ranges	EA
BLS38-47	OPTi ABV Refractometer % Mass w/w 0-35 and ABV (°Zeiss) (10-135) ranges	EA



Refractometers (Optical)

Portable, high quality, robust and easy to use zero adjust with lock. High precision clear scale. Easy clean prism. Low weight with comfortable hand grips and eyepiece.

Item	Description	Unit
BLS45-03	Eclipse Handheld optical refractometer. Sugar (Brix) 0-50%	EA
BLS45-07	Eclipse Handheld optical refractometer. Sugar (Brix) 0-32%	EA

Refractometers (Benchtop)

Decades of engineering excellence and experience have cemented Bellingham+Stanley refractometers as a favourite for reliable measurements in the drinks industry. The wide range enables the instrument for your budget and needs to be matched. Robust and reliable benchtop refractometers are complimented with light and strong digital handheld refractometers providing in-lab or in-situ measurement as required.

Item	Description	Unit
BLS19-00	Model RFM712-M 0-54 Brix 0.1 Brix resolution	EA
BLS19-10	Model RFM732-M 0-100 Brix 0.1 Brix resolution	EA
BLS19-20	Model RFM742-M 0-100 Brix 0.01 Resolution	EA



Refrigerators and Freezers - General Purpose Laboratory


A large range of general purpose refrigerators and freezers available from 151L compact underbench to large double door 1,400L models. Other sizes and model available.

Item Code	Description	Unit
TMO151R-AEW-TS	ES Series 151L +1c° to +10c° Underbench Refrigerator	EA
TMO288R-AEW-TS	ES Series 288L +1c° to +10c° Refrigerator	EA
TMO151F-AEW-TS	ES Series 151L -18c° to -25c° Underbench Freezer	EA
TMO232F-AEW-TS	ES Series 232L -18c° to -25c° Freezer	EA
TMOR700-GAEW-TS	GPS Series 700L +1c° to +11c° Glass Door Refrigerator	EA
TMOR14X-GAEW-TS	GPS Series 1400L +1c° to +11c° Glass Door Refrigerator	EA
TMOF700-SAEW-TS	GPS Series 700L -12c° to -25c° Solid Door Freezer	EA
TMOF14X-SAEW-TS	GPS Series 1400L -12c° to -25c° Solid Door Freezer	EA



Spectrophotometers (Cuvette)

Thermo Fisher Scientific offers a full range of cuvette spectrophotometers. Our systems are suitable for wine and grape juice chemical analysis and stability tests including glucose, malic acid, L-lactic acid and acetic acid monitoring.



Thermo Scientific Genesys 10S

The Thermo Scientific Genesys 10S spectrophotometers are available as visible light only or UV-Visible. They all come with a 6-position cell holder for rapid analysis, on-board control, USB method/data exporting, fast wavelength scanning and are simple to use.

Thermo Scientific Spectronic 200

The Thermo Scientific Spectronic 200 is a rugged, simple-to-use visible spectrophotometer for routine applications. It is suitable for basic wine and grape juice chemical analysis.

Item	Description	Unit
UNI840-209700	Genesys 10S UV-Vis: Includes 6 cell & 1 cell holders, scanning capability, 1.8nm bandwidth and Xenon flash lamp (Dual Beam, 190nm - 1100nm)	EA
UNI840-281600	Spectronic 200 Spectrophotometer with onboard software. Includes cell and test tube holder and 2 cuvette racks. Optional VISIONlite software for PC control	EA

Turbidimeters (Portable)

Orion AQUAfast IV Portable Advanced Turbidity Meter (AQ4500) has a dual source LED that allows readings that comply with EPA180.1 and ISO7027 methodologies for water. This dual light source has advantages in allowing measurements of both red and white wines. Supplied with carry case cells and batteries.

Item	Description	Unit
ORIAQ4500W	Orion AQUAfast IV portable advanced turbidity meter (AQ4500)	EA



Turbidimeters (Benchtop)

Easy operation with step-by-step menus to guide you through calibration, normal operation and diagnostics.

Ranges: 0 - 4000NTU (TL2300); 0 - 10000 NTU (TL2350). Resolution to three decimal places between 0 and 1.000 NTU.

Multiple unit selection: NTU or EBC.

Operator selectable ratio (for colour compensation) and non-ratio measurement modes. USB output featured for easy data transfer. Advanced ratio detectors for accurate low level measurement. Calibration using StablCal stabilised formazin turbidity standards.

Item	Description	Unit
HACLPV444.99.00210	TL2300 White light benchtop turbidimeter	EA
HACLPV444.99.00310	TL2350 White light benchtop turbidimeter	EA



Vortex

Low profile design with small footprint. Wide variety of optional accessories available to accommodate different vessel types.

Item	Description	Unit
BAR88880018	Vortex mixer with touch and continuous modes 0-3000 rpm analogue control	EA



Water Baths (Heated Circulated)

Sahara water baths have a choice of 5 thermostats and 8 stainless steel tanks, that can be configured to meet any temperature and volume requirements, giving you efficient heating for open and closed applications. Temperature range ambient+13°C to 200°C. Tank volume available 3-49L. Sahara water baths come with a 3 year guarantee.

Item	Description	Unit
HAK152-1211	SC100-S21 ambient+13°C to 100°C 19L tank with lid	EA
HAK152-1491	SC100-S49 ambient+13°C to 100°C 53L tank with lid	EA
HAK153-1211	SC150-S21 ambient+13°C to 150°C 19L tank with lid	EA



Water Baths (Refrigerated)

Arctic Refrigerated baths for applications that require low temperatures to -10°C.

Item	Description	Unit
HAK152-4101	Waterbath SC100-A10B -10 to +100°C 30L tank	EA
HAK152-4251	Waterbath SC100-A10B -25 to +100°C 21L tank	EA

Water Baths (Shaking)

Shaking water bath with maximum temperature 100°C and a horizontal motion.

Item	Description	Unit
BARTSSWB27	Precision shaking water bath 27L reciprocating digital, temperature range ambient+5°C up to 100°C. Interior 30 L x 69 W x 19 D cm exterior 36 L x 93 W x 24 D cm	EA



Water Baths (Unstirred)

Precision general purpose water baths are rugged, high performance baths that are designed to maintain water temperature from ambient to 100°C. Baths come with clear polycarbonate gable cover, diffuser tray and drain hose.

Item	Description	Unit
BARTSGP10	Precision model GP 10 unstirred water bath 10L digital display. Work area 30.1 L x 33.0 W x 15.0 H cm	EA
BARTSGP28	Precision model GP 28 unstirred water bath 28L digital display. Work area 29.7 L x 50.0 W x 20.0 H cm	EA



Water Purification Systems - Barnstead

Pure water is elemental to the success of your experiments. Our complete line of Barnstead water purification technologies includes solutions for critical and everyday application needs.

Item	Description	Unit
Type 1		
TSW50129872	Smart2Pure 3 UV Ultrapure Water Purification System - Type I water direct from tap water – suitable for HPLC / IC	EA
TSW50132373	MicroPure UV Ultrapure Water Purification System - Type I water from pretreated water – suitable for HPLC / IC	
TSW50131952	GenPure Pro UV Ultrapure Water Purification System - Type I Ultrapure	EA
Type 2		
TSW5013212	Pacific TII 3 Type II Water Purification System - Type II	EA
TSW50132123	Pacific TII 7 Type II Water Purification System - Type II	EA
TSW50132124	Pacific TII 12 Type II Water Purification System - Type II	EA
Type 3		
TSW50132385	Pacific RO 3 Water Purification System - Type III/RO	EA
TSW50132386	Pacific RO 7 Water Purification System - Type III/RO	EA
TSW50132387	Pacific RO 12 Water Purification System - Type III/RO	EA



Water Purification Systems - Elga

Elga provides water purification solutions for a variety of applications and industries. The PURELAB range provides innovation and flexibility to allow you to focus on your results, not your water.

Item	Description	Unit
Type 1		
ELGPC1GSCXM1	Chorus 1 General Science Water Purification System - Type I Ultrapure	EA
ELGLA755	Advanced Halo Dispenser	
Type 2		
ELGPC210DIXXM2	Chorus 2 10 L/hr Water Purification System - Type II (RO/DI)	EA
ELGPC220DIXXM2	Chorus 2 20 L/hr Water Purification System - Type II (RO/DI)	EA
Type 3		
ELGRO320XXM1	Chorus 3 20 L/hr Water Purification System - Type III/RO	EA
ELGRO330XXM1	Chorus 3 30 L/hr Water Purification System - Type III/RO	EA



Filtration

Membrane Filters	43
Filter Holders	43
Glass Fibre Filter Paper	43
Qualitative Filter Paper	44
Syringe Filters	45

A critical first step in so many wine analyses, our ranges include Whatman, Fisherbrand and Millipore brands of filter papers, glass fibre, membranes and syringe filters. Apparatus are also available, from the smallest to largest sample volumes.

Membrane Filters

Millipore and Whatman membrane filters for microbiological testing. Sterilised & individually packaged in peel back envelopes.

Item	Description	Unit
MILHAWG04700	Membrane filter 0.45µm 47mm grid MCE	100pcs
MILHAWG047S6	Membrane filter 0.45µm 47mm S-Pak MCE	600pcs
MILHAWP04700	Membrane filter 0.45µm 47mm plain MCE	100pcs
MILHVLPO4700	Membrane filter durapore 0.45µm 47mm plain PVDF	100pcs
WHA7141-154	Membrane filter 0.45µm 47mm MCE gridded sterile bulk pack	1000pcs
Cellulose Nitrate - Non Sterile		
SAT1132725N	Filter membrane CN 0.2µm 25mm	100pcs
SAT1132747N	Filter membrane CN 0.2µm 47mm	100pcs
SAT1130625N	Filter membrane CN 0.45µm 25mm	100pcs
SAT1130647N	Filter membrane CN 0.45µm 47mm	100pcs
Cellulose Nitrate - Sterile Individually Packed		
SAT1130625ACN	Filter membrane CN 0.45um 25mm	100pcs
SAT1130647ACN	Filter membrane CN 0.45µm 47mm	100pcs
White membrane with black grid - Sterile Individually Packed		
SAT1140747ACN	Filter membrane CN 0.2um 47mm	100pcs
SAT1140647ACN	Filter membrane CN 0.45um 47mm	100pcs
SAT1140747ACR	Filter membrane CN 0.2um 47mm	1000pcs
SAT1140647ACR	Filter membrane CN 0.45um 47mm	1000pcs

Filter Holders

Item	Description	Unit
MILXX1004700	Filter holder glass 47mm, 300mL capacity, with frit glass filter, includes clamp and stopper	EA
MILXX6200006P	Filter forceps blunt end stainless steel	3pcs
SCH2120154	Filter flask borosilicate 1 litre with glass side arm connection	EA

Glass Fibre Filter Paper

Item	Description	Unit
Fisherbrand G8 – Equivalent to Whatman GF/D		
FSB09-804-21D	Filter paper glass fibre G8 2.1cm 2.5µm	100pcs
FSB09-804-25D	Filter paper glass fibre G8 2.5cm 2.5µm	100pcs
FSB09-804-42D	Filter paper glass fibre G8 4.25cm 2.5µm	100pcs
FSB09-804-55D	Filter paper glass fibre G8 5.5cm 2.5µm	100pcs
FSB09-804-70D	Filter paper glass fibre G8 7cm 2.5µm	100pcs
Fisherbrand G6 – Equivalent to Whatman GF/A		
FSB09-804-21A	Filter paper glass fibre G6 2.1cm 1.6µm	100pcs
FSB09-804-24A	Filter paper glass fibre G6 2.4cm 1.6µm	100pcs
FSB09-804-42A	Filter paper glass fibre G6 4.25cm 1.6µm	100pcs
FSB09-804-55A	Filter paper glass fibre G6 5.5cm 1.6µm	100pcs



Glass Fibre Filter Paper (Continued)

Item	Description	Unit
FSB09-804-70A	Filter paper glass fibre G6 7cm 1.6µm	100pcs
FSB09-804-90A	Filter paper glass fibre G6 9cm 1.6µm	100pcs
FSB09-804-110A	Filter paper glass fibre G6 11.0cm 1.6µm	100pcs
FSB09-804-125A	Filter paper glass fibre G6 12.5cm 1.6µm	100pcs
FSB09-804-150A	Filter paper glass fibre G6 15cm 1.6µm	100pcs
WHA1820-047	Filter glass microfibre GF/A 1.6µm 4.7cm	100pcs
WHA1820-110	Filter glass microfibre GF/A 1.6µm 11cm	100pcs

Fisherbrand G4 – Equivalent to Whatman GF/C

FSB09-804-21C	Filter paper glass fibre G4 2.1cm 1.2µm	100pcs
FSB09-804-24C	Filter paper glass fibre G4 2.4cm 1.2µm	100pcs
FSB09-804-42C	Filter paper glass fibre G4 4.25cm 1.2µm	100pcs
FSB09-804-55C	Filter paper glass fibre G4 5.5cm 1.2µm	100pcs
FSB09-804-70C	Filter paper glass fibre G4 7cm 1.2µm	100pcs
FSB09-804-90C	Filter paper glass fibre G4 9cm 1.2µm	100pcs
WHA1822-025	Filter glass microfibre GF/C 1.2µm 2.5cm	100pcs
WHA1822-047	Filter glass microfibre GF/C 1.2µm 4.7cm	100pcs
WHA1822-150	Filter glass microfibre GF/C 1.2µm 15cm	100pcs

Fisherbrand G2 – Equivalent to Whatman GF/B

FSB09-804-21B	Filter paper glass fibre G2 2.1cm 1.0µm	100pcs
FSB09-804-55B	Filter paper glass fibre G2 5.5cm 1.0µm	100pcs
WHA1821-025	Filter glass microfibre GF/B 1.0µm 2.5cm	100pcs



Qualitative Filter Paper

Item	Description	Unit
WHA1001-047	Filter paper qualitative no.1 11µm 4.7cm	100pcs
WHA1001-070	Filter paper qualitative no.1 11µm 7cm	100pcs
WHA1001-090	Filter paper qualitative no.1 11µm 9cm	100pcs
WHA1001-110	Filter paper qualitative no.1 11µm 11cm	100pcs
WHA1001-125	Filter paper qualitative no.1 11µm 12.5cm	100pcs
WHA1001-150	Filter paper qualitative no.1 11µm 15cm	100pcs

Fisherbrand P8 – Equivalent to Whatman No.4

FSB09-795-AA	Filter paper P8 2.5cm 20-25µm	100pcs
FSB09-795-A	Filter paper P8 5.5cm 20-25µm	100pcs
FSB09-795-B	Filter paper P8 7cm 20-25µm	100pcs
FSB09-795-C	Filter paper P8 9cm 20-25µm	100pcs
FSB09-795-D	Filter paper P8 11cm 20-25µm	100pcs
FSB09-795-E	Filter paper P8 12.5cm 20-25µm	100pcs
FSB09-795-F	Filter paper P8 15cm 20-25µm	100pcs
FSB09-795-G	Filter paper P8 18.5cm 20-25µm	100pcs
FSB09-804	Filter paper P8 19cm 20-25µm	100pcs
FSB09-795-H	Filter paper P8 20.5cm 20-25µm	100pcs



Qualitative Filter Paper (Continued)

Item	Description	Unit
FSB09-795-J	Filter paper P8 24cm 20-25µm	100pcs
FSB09-795-K	Filter paper P8 33cm 20-25µm	100pcs
FSB09-800	Filter paper P8 50x50cm 20-25µm	500pcs
FSB09-802-1B	Filter paper P8 46x57cm 20-25µm	100pcs
WHA1004-047	Filter paper qualitative no. 4 20-25µm 4.7cm	100pcs
WHA1004-125	Filter paper qualitative no. 4 20-25µm 12.5cm	100pcs
Fisherbrand P2 – Equivalent to Whatman No.5		
FSB09-803-5A	Filter paper P2 5.5cm 1-5µm	100pcs
FSB09-803-5B	Filter paper wP2 7cm 1-5µm	100pcs
FSB09-803-5C	Filter paper P2 9cm 1-5µm	100pcs
FSB09-803-5E	Filter paper P2 12.5cm 1-5µm	100pcs
FSB09-803-5F	Filter paper P2 15cm 1-5µm	100pcs
WHA1005-110	Filter paper qualitative no.5 2.5µm 11cm	100pcs
WHA1005-150	Filter paper qualitative no.5 2.5µm 15cm	100pcs



Syringe Filters

Syringe filters offer a reliable and fast means of clarifying and preparing sample volumes of 1mL-200mL. A wide variety of membrane types are available to suit your application.

Item	Description	Unit
FSBSE4M081I99	Syringe filter 25mm MCE 0.45µm non sterile	100pcs
FSBSE4M018I09	Syringe Filter 33mm Nylon 0.2µm non sterile	50pcs
FSBSE4M019I09	Syringe Filter 33mm Nylon 0.2µm non sterile	150pcs
FSBSE4J020I09	Syringe Filter 33mm Nylon 0.45µm non sterile	50pcs
FSBSE4J021I09	Syringe Filter 33mm Nylon 0.45µm non sterile	150pcs
FSBSE4M036I09	Syringe Filter 33mm Nylon 0.45µm non sterile	1500pcs
FSBSE2M787H5	Syringe Filter 25mm PTFE 0.2µm non sterile	50pcs
FSBSE4M039I00	Syringe Filter 25mm PTFE 0.2µm non sterile	150pcs
FSBSE2M780H5	Syringe Filter 25mm PTFE 0.45µm non sterile	50pcs
FSBSE4M040I00	Syringe Filter 25mm PTFE 0.45µm non sterile	150pcs
THCCH4525-PES	Syringe Filter 25mm 0.45µm PES	100pcs
THCCH4525-PVH	Syringe Filter 25mm 0.45µm PVDF hydrophobic	100pcs
THCCH4525-CA	Syringe Filter 25mm 0.45µm CA	100pcs
THCCH4525-NN	Syringe Filter 25mm 0.45µm Nylon	100pcs
THCCH4525-NP	Syringe Filter 25mm 0.45µm PTFE	100pcs
THCCH4525-PP	Syringe Filter 25mm 0.45µm Polypropylene	100pcs
THCCH4525-PV	Syringe Filter 25mm 0.45µm PVDF hydrophilic	100pcs
SAT16534K	Syringe filter Minisart NML 0.2µm SFCA 28mm	50pcs
SAT16555K	Syringe filter Minisart NML 0.45µm SFCA 28mm	50pcs
SAT16534Q	Syringe filter Minisart NML 0.2µm SFCA 28mm	500pcs
SAT16555Q	Syringe filter Minisart NML 0.45µm SFCA 28mm	500pcs



Find Syringes on page 25

Microbiology & Hygiene Monitoring

Dehydrated Culture Media	47
Gas Generating Systems	47
Hygiene Monitoring	47
Prepared Media	48
Swabs	48
Test Kits (Allergen)	49

All your needs covered with our extensive range of microbiology products, including hygiene monitoring, pre-prepared media and enzyme test kits.

Dehydrated Culture Media

Oxoid dehydrated culture media for the determination of the microbiological flora in the fermentation processes.

Item	Description	Unit
CM0309B	WL Nutrient Agar	500g
CM0501B	WL Nutrient Broth	500g
CM0361B	MRS Agar	500g
SR0222C	Cyclohexamide Solution 0.1%	100mL



Gas Generating Systems

Oxoid gas generating systems for creation of the optimal growth conditions for anaerobes, microaerophics and CO₂ dependent organisms.

Item	Description	Unit
AG0025A	AnaeroJar 2.5L (fully assembled) holds up to 12x 90mm plates	EA
AN0025A	AnaeroGen 2.5L Anaerobic atmosphere generation sachets	10pcs
AN0035A	AnaeroGen 3.5L Anaerobic atmosphere generation sachets	10pcs
AN0010C	AnaeroGen compact complete with 10 sachets and pouches	10pcs
BR0055B	Anaerobic indicator	100pcs



Hygiene Monitoring

Neogen AccuPoint® Advanced: Hygiene Monitoring System is a lightweight, handheld system that validates the effectiveness of a sanitation program by detecting ATP from food residues and microorganisms present on surfaces and in liquids. The system utilises easy-to-use AccuPoint samplers, which automatically mix prepared reagents to quickly generate results.

Neogen AccuClean® is intended to be used as a simple method to determine the cleanliness of food contact surfaces and equipment. After 10 minutes, the sampling tip is compared to a color chart to gauge the level of protein detected.

Item	Description	Unit
NEO9902	AccuPoint Advanced Hygiene Monitoring System	EA
NEO9905	AccuPoint Advanced Surface Sampler	100pcs
NEO9906	AccuPoint Advanced Water Sampler	100pcs



Microbiology & Hygiene Monitoring


Hygiene Monitoring (Continued)

Item	Description	Unit
NEO9907	AccuPoint Advanced Access Sampler	100pcs
NEO9670	AccuClean Sanitation Test Sampler	100pcs
DS0147A	Dip Slide TTC Red Spot	10pcs
DS0155A	Dip Slide TTC Red Spot / Malt Extract Agar	10pcs
EM0100S	Oxoid Air Sampler with Stainless Steel Head for 90mm plates	EA

Prepared Media

Oxoid prepared plated and tube media are manufactured in our NATA accredited facility. 3M™ Petrifilm™ AQUA Plates are sample-ready plates that replace conventional agar, petri dishes, media pads and disposable filter funnels used in the microbial testing of water. Colitag rapid presence/absence test for *E.coli* & coliform in water.

Item	Description	Unit
PP2145	Plate Count Agar	10pcs
PP2159	R2A Agar	10pcs
PP2311	MRS Agar	10pcs
PP2358	Lysine Medium	10pcs
PP2670	Chromogenic coliform agar (ISO)	10pcs
PP2296	WLNA + Cycloheximide	10pcs
PP2136	WLDA	10pcs
PP2137	WLNA	10pcs
TM1864	Brettanomyces/Dekkara Broth 15mL	75pcs
3M6450	3M™ Petrifilm™ AQUA Heterotrophic count plate	100pcs
3M6457	3M™ Petrifilm™ AQUA Coliform count plate	50pcs
3M6408	3M™ Petrifilm™ AQUA Yeast & Mould count plate	100pcs
3M6418	3M™ Petrifilm™ AQUA Enterobacteriaceae plate	50pcs
CPI4600-0011	Colitag test kit	20pcs
CPI4600-0013	Colitag test kit	100pcs
CPI4600-0121C	Colitag 120mL sample vessels with sodium thiosulphate	100pcs
CPI4600-FILM5	Sealing Film	5pcs
CPI4600-ST	MPN Plate Sealing Tool Roller	EA
CPI46001	MPN Plate	100pcs



Swabs

Item Code	Description	Unit
LABMC150NS	Swab non sterile 150mm cotton tip	100pcs
TM4222	Environmental Flexi Swab	100pcs
MWE102-800	Swab sterile dry plastic shaft	800pcs

Test Kits (Allergen)

Elisa Systems provide allergen test kits to meet the new EU allergens labelling requirements.

Item	Description	Unit
ESEGG-48	Elisa Systems Enhanced Egg Allergen Kit	48pcs
ESCASPRD	Elisa Systems Casein Allergen Kit	48pcs
ESADDSOL	Elisa Systems Allergen Kit Extraction Solution for Polyphenols in red wine samples 20X concentrate to make a final volume of 500mL	EA



Wine analysis resources cellar

Understanding the demands of wine safety testing can be a complex task for the laboratory with numerous analytes and residues to monitor. Our separation and detection technologies, combined with experienced applications competence, and our best suited chemistries provide ideal solutions for the analysis of wine. Each library is a targeted collection of scientific applications, videos, webinars, and resources based on your interests.

Visit our wine analysis resources cellar at [thermofisher.com/winetesting](https://www.thermofisher.com/winetesting)

Reagent Kits & Analysers

Instrument: Automated Photometric Analyser	51
Instrument: Automatic Analyser (FTIR-UV-VIS)	51
Instrument: GC and GC/MS	52
Instrument: Ion Chromatography	52
Instrument: Isotope Ratio Mass Spectrometry	52
Instrument: Liquid Chromatography	53
Instrument: Trace elemental	53
Pathogen Detection	55
Pesticide Analysis	56
System Reagents	58
System Consumables	59

Increasing productivity by automating testing and quality monitoring.

Automated Photometric Analyser

Thermo Scientific Gallery and Gallery Plus discrete analysers are easy to use, automated systems that allow laboratories to simplify their testing with the dual benefits of time and cost savings. All necessary analysis steps are automated, providing true walk-away time for the operator. Applications include food, beverage, water, soil, and industrial quality control.



In wine production, sugars, acids, alcohol, and sulfite analysis are important parameters for production and quality control. Gallery analysers require no advanced preparation time and samples can be added at any time without interrupting the current run.

Both Gallery instruments provide an integrated platform for two measurement techniques, photometric and electrochemical (ECM) which can be run simultaneously.

Item Code	Description	Unit
CDK98610001-2	Gallery with wine filter set	EA
CDK98611001-2	Gallery with ECM (pH electrode) and Wine Filter set	EA
CDK98620001-2	Gallery Plus with Wine Filter set	EA
CDK98621001-2	Gallery Plus with ECM (pH electrode) with wine filter set	EA

Automatic Analyser (FTIR - UV - VIS)

The Bacchus 3 Automatic Analyser (FTIR - UV - VIS) is the only analyser on the marketplace that provides simultaneously the integral IR, UV and visible absorbance spectra; it is therefore particularly well adapted to cover most wine analysis situations, including; Maturity control at reception bay level; Must and fermenting must; Finished wines (dry, sweet, sparkling, fortified wines).



Item Code	Description	Unit
BACCHUS3	Bacchus3 - Base	EA
BACCHUS3-AUTO	Bacchus3 – Base with Autosampler	EA
BACCHUS3-MULTISPEC	Bacchus3 - UV/Vis spectrometer module (optional)	EA

Reagent Kits & Analysers

Gas Chromatography (GC) and Gas Chromatography–Mass Spectrometry (GC-MS) Systems

Taints analysis

Complete solutions for the analyses of beverages including contaminants, flavors and FAMES. Advanced sampling capabilities, the unique instant-connect injectors and detectors, and no-vent MS column replacement provide maximum flexibility and increased productivity by minimizing maintenance downtime.

Item Code	Description	Unit
Contact us	Thermo Scientific TRACE 1310 GC + Thermo Scientific TriPlus 300 Headspace AutoSampler	EA
Contact us	Thermo Scientific ISQ Series Quadrupole GC-MS HPIC	EA
Contact us	Thermo Scientific TSQ Triple Quadrupole GC-MS/MS	EA
Contact us	Thermo Scientific Q Exactive Mass Spectrometer	EA
Contact us	Thermo Scientific Exactive GC Orbitrap GC-MS/MS	EA



Ion Chromatography

Determination of High Sugar Concentrations

Quantify levels of glucose and fructose present in wine using the Thermo Scientific™ Dionex™ Integrion™ HPIC™ system.

Organic Acids by High-Pressure IC

For more demanding beverage applications, the best solution is the Thermo Scientific™ Dionex™ ICS-5000+ Reagent-Free™ HPIC™ system which can be used to resolve a large number of organic acids and inorganic anions using hydroxide gradient elution.

Item Code	Description	Unit
Contact us	Thermo Scientific Dionex Integrion HPIC	EA
Contact us	Thermo Scientific Dionex ICS-5000+ Reagent-Free HPIC	EA



Isotope Ratio Mass Spectrometry

Isotope Analysis of Wine

The measurement of the $^{18}\text{O}/^{16}\text{O}$ ratio of water is one of the most important applications of isotope ratio mass spectrometry. Authenticity control of wine demands an analytical solution that can provide high precision, high accuracy and high throughput analysis. ^{18}O isotopic analysis of wine has been a topic of increasing interest, and is now the subject of EU regulation and other countries.

Item Code	Description	Unit
Contact us	Thermo Scientific GasBench II Isotope Ratio Mass Spectrometer (IRMS)	EA



Reagent Kits & Analysers

Liquid Chromatography

Analysis and quantification of polyphenolics

Polyphenols such as catechins and phenolic acids are aromatic compounds found in all fruits and vegetables. Red grape skins are particularly rich in polyphenols which makes red wine a great source of these molecules. But to determine the quality and taste of the wine, the levels of polyphenols must be quantified.

Pesticide residue analysis

Monitoring pesticides residue is incredibly important because they are generally toxic and damage human health with long-term exposure.

Both these applications can be determined using the Thermo Scientific™ Vanquish™ UHPLC System which leverages the latest, patented LC technology to boost uptime, extend robustness and increase result reliability.

UHPLC with Mass Spectrometry

Red wine is a highly complex mixture of molecules including polyphenols, carbohydrates, proteins and flavanols, reported to have antioxidant properties. Ultra High Pressure Liquid Chromatography (UHPLC) coupled to the Thermo Scientific™ Q Exactive™ Focus Hybrid Quadrupole-Orbitrap™ mass spectrometer permits simultaneous detection and quantitation of these particularly important wine constituents.

Item Code	Description	Unit
Contact us	Thermo Scientific Vanquish UHPLC System	EA
Contact us	Thermo Scientific Q Exactive Focus Hybrid Quadrupole-Orbitrap mass spectrometer	EA



Trace elemental

Analysis of heavy metals

Heavy metals such as copper, zinc and iron are commonly found in wines. They are problematic and can dramatically affect sensory and shelf-life quality if present at too high a level. There are many contributing factors that determine the metal content in wine including; soil and type of vineyard to name a few. Heavy metal testing before bottling is always recommended. The prescribed method of analysis is Flame Atomic Absorption Spectroscopy (FAAS) and the Thermo Scientific™ iCE™ 3400 Atomic Absorption Spectrometer suited for this role.

Analysis of trace elements

Further, you can also use the Thermo Scientific™ iCAP™ RQ ICP-MS and the Thermo Scientific™ iCAP™ 7400 ICP-OES to analyse certain parameters to determine the quality of a wine such as key toxic and non-toxic elements, as well as, nutritional components.



Item Code	Description	Unit
Contact us	Thermo Scientific iCE 3400 Atomic Absorption Spectrometer	EA
Contact us	Thermo Scientific iCAP RQ ICP-MS	EA
Contact us	Thermo Scientific iCAP 7400 ICP-OES	EA



What vintage are your analytical tools?

At Thermo Fisher Scientific, we won't let equipment financing stand between you and your next great discovery. Thermo Fisher is your one-stop partner for the best laboratory products and analytical technologies available, plus the unique financing options you need to accelerate success in science or industry. These innovative options can help meet your company's budgetary needs and bottom-line goals.

Find out more at thermofisher.com.au/finance

ThermoFisher
SCIENTIFIC

Reagent Kits & Analysers

Pathogen Detection

Bioreba was the first company worldwide to produce and commercialise ELISA reagents for plant virus diagnostics. Choose from a range of high quality products to address the detection of a wide range of viruses that effect vineyards and their quality of vines. In addition, choose from a range of extraction bags available in four different sizes, with or without a synthetic intermediate layer.



Item Code	Description	Unit
BRB120372	Arabis mosaic virus (ArMV) 960 Assay Complete Kit	EA
BRB120375	Arabis mosaic virus (ArMV) 480 Complete Kit	EA
BRB120272	Arabis mosaic virus + Grapevine fanleaf virus (ArMV+GFLV) 960 Complete Kit	EA
BRB120275	Arabis mosaic virus + Grapevine fanleaf virus (ArMV+GFLV) 480 Complete Kit	EA
BRB120672	Grapevine fleck virus (GFkV) 960 Complete Kit	EA
BRB120675	Grapevine fleck virus (GFkV) 480 Assays Complete Kit	EA
BRB120472	Grapevine fanleaf virus (GFLV) 960 Assays Complete Kit	EA
BRB120475	Grapevine fanleaf virus (GFLV) 480 Assays Complete Kit	EA
BRB120172	Grapevine leafroll-assoc. virus 1 (GLRaV-1) 960 Assays Complete Kit	EA
BRB120175	Grapevine leafroll-assoc. virus 1 (GLRaV-1) 480 Assays Complete Kit	EA
BRB120972	Grapevine leafroll-assoc. virus 1+3 (GLRaV-1+3) 960 Assays Complete Kit	EA
BRB120975	Grapevine leafroll-assoc. virus 1+3 (GLRaV-1+3) 480 Assays Complete Kit	EA
BRB120872	Grapevine leafroll-assoc. virus 2 (GLRaV-2) 960 Assays Complete Kit	EA
BRB120875	Grapevine leafroll-assoc. virus 2 (GLRaV-2) 480 Assays Complete Kit	EA
BRB122372	Grapevine leafroll-assoc. virus 3 (GLRaV-3) 960 Assays Complete Kit	EA
BRB122375	Grapevine leafroll-assoc. virus 3 (GLRaV-3) 480 Assays Complete Kit	EA
BRB122572	Grapevine leafroll generic 4-9 (GLRaV-4-9) 960 Assays Complete Kit	EA
BRB122575	Grapevine leafroll generic 4-9 (GLRaV-4-9) 480 Assays Complete Kit	EA
BRB122472	Grapevine leafroll-assoc. virus 6 (GLRaV-6) 960 Assays Complete Kit	EA
BRB122475	Grapevine leafroll-assoc. virus 6 (GLRaV-6) 480 Assays Complete Kit	EA
BRB121072	Grapevine leafroll-assoc. virus 7 (GLRaV-7) 960 Assays Complete Kit	EA
BRB121075	Grapevine leafroll-assoc. virus 7 (GLRaV-7) 480 Assays Complete Kit	EA
BRB122272	Grapevine virus A (GVA) 960 Assays Complete Kit	EA
BRB122275	Grapevine virus A (GVA) 480 Assays Complete Kit	EA



Extraction Bags

Item Code	Description	Unit
BRB430100A	Extraction bags universal 12 x 15 cm with synthetic intermediate layer for all kinds of samples	EA
BRB480100A	Extraction bags universal U-form 12 x 12 cm with synthetic layer for testing seeds	EA
BRB420100A	Extraction bags universal long 15 x 50cm without intermediate layer for testing long leaves	EA
BRB470100A	Extraction bags universal long 15 x 28cm heavy duty with synthetic intermediate layer	EA



Reagent Kits & Analysers

Pesticide Analysis - QuEChERS Dispersive SPE Extraction and Clean-up

QuEChERS (Pronounced 'catchers' - Quick Easy Cheap Effective Rugged and Safe) is a dispersive Solid Phase Extraction (SPE) technique for extracting multi-residues of pesticides from fruits and vegetables. The QuEChERS method is a two-step process: extraction followed by clean up. Pre-packed extraction kits contain pre-weighed salts and are supplied in 50mL PP centrifuge tubes. Sample clean-up kits are supplied in a variety of tube formats and bed weights to accommodate large and small sample sizes.



Item Code	Description	Unit
QuEChERS Extraction Products		
THC60105-210	50mL Tube PP with 6g Anhydrous Magnesium Sulphate & 1.5g Anhydrous Sodium Acetate 250pcs	EA
THC60105-211	50mL tube with 4g Anhydrous Magnesium Sulphate & 1g Anhydrous Sodium Chloride 250pcs	EA
THC60105-212	50mL tube with 6g Anhydrous Magnesium Sulphate 1.5g Sodium Citrate Dihydrate & 0.75g Disodium Citrate Sesquihydrate 250pcs	EA
THC60105-216	50mL tube with 1g Anhydrous Magnesium Sulphate 1g Sodium Citrate Dihydrate & 0.5g Disodium Citrate Sesquihydrate 250pcs	EA
QuEChERS Clean-up Products		
THC60105-219	2mL Tube with 150mg Anhydrous Magnesium Sulphate 25mg PSA 100pcs	EA
THC60105-220	2mL Tube with 150mg Anhydrous Magnesium Sulphate 25mg PSA & 25mg C18 100pcs	EA
THC60105-221	2mL Tube with 150mg Anhydrous Magnesium Sulphate 25mg PSA & 2.5mg Graphitized Carbon 100pcs	EA
THC60105-222	2mL Tube with 150mg Anhydrous Magnesium Sulphate 25mg PSA & 7.5mg Graphitized Carbon 100pcs	EA
THC60105-202	2mL Tube with 150mg Anhydrous Magnesium Sulphate 50mg PSA & 50mg Graphitized Carbon 100pcs	EA
THC60105-203	2mL Tube with 150mg Anhydrous Magnesium Sulphate 50mg PSA 100pcs	EA
THC60105-204	2mL Tube with 150mg Anhydrous Magnesium Sulphate 50mg PSA & 50mg Endcapped C18 100pcs	EA
THC60105-223	2mL Tube with 150mg Anhydrous Magnesium Sulphate 50mg PSA & 50mg Graphitized Carbon 50mg 100pcs	EA
THC60105-206	15mL Tube with 900mg Anhydrous Magnesium Sulphate 300mg PSA & 150mg C18 50pcs	EA

Reagent Kits & Analysers

Item Code	Description	Unit
QuEChERS Clean-up Products (Continued)		
THC60105-205	15mL Tube with 900mg Anhydrous Magnesium Sulphate 300mg PSA & 150mg Graphitized Carbon 50pcs	EA
THC60105-224	15mL Tube with 1200mg Anhydrous Magnesium Sulphate 400mg PSA 50pcs	EA
THC60105-225	15mL Tube with 1200mg Anhydrous Magnesium Sulphate 400mg PSA & 400mg Endcapped C18 50pcs	EA
THC60105-226	15mL Tube with 1200mg Anhydrous Magnesium Sulphate 400mg PSA 400mg Endcapped C18 & 400mg Graphitized Carbon 50pcs	EA
THC60105-227	15mL Tube with 900mg Anhydrous Magnesium Sulphate 150mg PSA & 150mg C18 50pcs	EA
THC60105-230	15mL Tube with 150mg Anhydrous Magnesium Sulphate 50mg PSA & 50mg Graphitized Carbon 50pcs	EA
THC60105-217	15mL Tube with 900mg Anhydrous Magnesium Sulphate 150mg PSA & 45mg Graphitized Carbon 50pcs	EA
THC60105-218	15mL Tube with 900mg Anhydrous Magnesium Sulphate 150mg PSA & 15mg Carbon 50pcs	EA
THC60105-215	15mL Tube with 900mg Anhydrous Magnesium Sulphate & 150mg PSA 50pcs	EA
THC60105-231	15mL Tube with 150mg Anhydrous Magnesium Sulphate 300mg PSA & 150mg Chloro Filter 50pcs	EA
THC60105-213	15mL Tube with 750mg Anhydrous Magnesium Sulphate 250mg PSA 250mg Endcapped C18 & 250mg Graphitized Carbon 50pcs	EA
THC60105-214	15mL Tube with 900mg Anhydrous Magnesium Sulphate & 300mg PSA 50pcs	EA
QuEChERS Multipacks		
THC60105-332	4000mg Magnesium Sulphate 1000mg Sodium Chloride	EA
THC60105-333	4000mg Magnesium Sulphate 1000mg Sodium Chloride 500mg Sodium Citrate dibasic sesquihydrate 1000mg Sodium Citrate tribas	EA
THC60105-334	4000mg Magnesium Sulphate 1g Sodium Acetate	EA
THC60105-335	6000mg Magnesium Sulphate 1500mg Sodium Acetate	EA
THC60105-336	6000mg Magnesium Sulphate 1500mg Sodium Chloride	EA
THC60105-337	6000mg Magnesium Sulphate 1500mg Magnesium Sulphate 1500mg Sodium citrate dihydrate 750mg disodium citrate sesquihydrate	EA
THC60105-338	8000mg Magnesium Sulphate 200mg of Sodium Chloride	EA
THC60105-339	8000mg Magnesium Sulphate 3500mg of Sodium Chloride	EA



Reagent Kits & Analysers

Pesticide Analysis (Continued)

Item Code	Description	Unit
QuEChERS Mylar Pouches		
THC60105-340	Mylar pouch: 4000mg Magnesium Sulphate 1000mg Sodium Chloride - (no tubes)	EA
THC60105-341	Mylar pouch: 6000mg Magnesium Sulphate 1500mg CH ₃ COONa (anhydrous) - (no tubes)	EA
THC60105-342	Mylar pouch: 6000mg Magnesium Sulphate 1500mg Sodium Chloride - (no tubes)	EA
THC60105-343	Mylar pouch: 6000mg Magnesium Sulphate 1500mg Sodium Chloride 1500mg Na ₃ Cit2H ₂ O 0.75g Na ₂ HCit - (no tubes)	EA
THC60105-344	Mylar pouch 4000mg Magnesium Sulphate 1000mg Sodium Chloride 500mg Na ₂ Cit - (no tubes)	EA
THC60105-345	1.2g Sodium Chloride in a bag	EA
THC60105-346	4g disodium Sulphate and 0.5g Magnesium Sulphate	EA

System Reagents (for Gallery & Konelab Arena systems)

Item	Description	Unit
CDK984347	Acetaldehyde 250 tests	EA
CDK984303	Acetic Acid 1000 tests	EA
CDK984318	Acetic Acid 300 tests	EA
CDK984342	Alpha-Amino Nitrogen (NOPA) 300 tests	EA
CDK984320	Ammonia 300 tests	EA
CDK984635	L-Ascorbic Acid 250 tests	EA
CDK984361	Calcium 350 tests	EA
CDK984327	Citric Acid 250 tests	EA
CDK984628	Copper 300 tests	EA
CDK984302	D-Fructose 300 tests	EA
CDK984710	D-Gluconic Acid 250 tests	EA
CDK984304	D-Glucose 300 tests	EA
CDK984313	D-Glucose + D-Fructose 1000 tests	EA
CDK984314	D-Glucose + D-Fructose 300 tests	EA
CDK984301	D-Glucose + D-Fructose + Sucrose 1000 tests	EA
CDK984317	D-Glucose + D-Fructose + Sucrose 300 tests	EA



Reagent Kits & Analysers

System Reagents (for Gallery & Konelab Arena systems) (Continued)

Item	Description	Unit
CDK984300	Ethanol 350 tests	EA
CDK984316	Glycerol 300 tests	EA
CDK984308	L-Lactic Acid 300 tests	EA
CDK984310	L-Malic Acid 300 tests	EA
CDK984311	L-Malic Acid 1000 tests	EA
CDK984360	Magnesium 350 tests	EA
CDK984348	Oxalic Acid 250 tests	EA
CDK984349	pH (Colorimetric) (330)	EA
CDK984634	SO ₂ Free 430 tests	EA
CDK984345	SO ₂ Total 300 tests	EA
CDK984312	Sucrose (Total Glucose) 300 tests	EA
CDK984632	Total Acids (Wine pH7) 400 tests	EA
CDK984326	Total Iron 850 tests	EA
CDK984346	Total Polyphenol 280 tests	EA
CDK984328	Total Protein (Biuret)	EA
CDK984321	Urea (Ammonia) 775 tests	EA

System Consumables

Item	Description	Unit
CDK986540	DECACELL cuvettes, 5400 tests (for Gallery)	EA
CDK984000	Multicell cuvettes, 40 x 25pcs (for Konelab Arena)	EA
CDK989221	2mL sample cups	1000pcs
BPL8365	Blood/Culture Caps 10 12 13 16 18mm Green Hexa-Flex Safety Tube Caps	500pcs
CDK984619	7mL sample tubes	1000pcs
CDK984030	Washing solution 4.5%	4 x 20mL
CDK981456	20mL reagent bottles	16pcs
LBS40001N	LabServe Cup 2.5mL for Technicon AA clear	5000pcs



Cellular & Molecular Tools for Wine Analysis

Cell Counter	61
Cell Imaging Systems	62
Flow Cytometer	62
Nucleic Acid Purification	63
Real Time PCR Instruments & Assays	63

For products in this section please place your order with the following legal entity:

Life Technologies Australia Pty Ltd

Life Technologies is now part of Thermo Fisher Scientific. Although we share the same Company name we are still operating as two separate legal entities.

Please ensure all products in this section are ordered through our legal entity: Life Technologies Australia Pty Ltd.

There is no change to how you currently do business with either legal entity. Our service and support teams remain the same and are focused in providing specific assistance for those products sold by each entity. If you have any questions or require further detail, please contact our customer care team.

To order Life Science Solutions products, in this section please place your order with:

Life Technologies Australia Pty Ltd.

Customer Service: 1800-636-327

Technical Support: 1800-636-327

Fax: 1800-143-363

Email (Orders): orders.anz@lifetech.com

Email (Inquires): mail.anz@lifetech.com

Sales Representatives: thermofisher.com/salesrep

For all other products in this catalogue please place your order with:

Thermo Fisher Scientific Australia Pty Ltd.

Customer Service: 1300-735-292

Technical Support: 1300-736-767

Fax: 1800-067-639

Email (Orders): ordersau@thermofisher.com

Email (Inquires): auinfo@thermofisher.com

Website: thermofisher.com.au

Cell Counter

With the option for a reusable slide and fluorescence capabilities—brightfield and two user-changeable fluorescence channels—the Invitrogen™ Countess™ II FL Automated Cell Counter can count cells as small as 7µm, monitor fluorescent protein expression, and measure cell viability.

Find out more at thermofisher.com/countess

Item	Description	Unit
AMQAF1000	Countess® II FL Automated Cell Counter	EA
AMQAX1000	Countess® II Automated Cell Counter	EA

Only available for order through legal entity Life Technologies Australia Pty Ltd



Cellular & Molecular Tools for Wine Analysis

Cell Imaging Systems

Designed to eliminate the complexities of microscopy without compromising performance, the EVOS® line of cell imaging systems makes cell imaging accessible to almost every lab and budget.

Find out more at thermofisher.com/evos

Item	Description	Unit
AMEX1000	EVOS® XL Core Imaging System	EA
AME3300	EVOS® XL Imaging System	EA
4471136	FLoid® Cell Imaging Station	EA
AMF4300	EVOS® FL Imaging System	EA
AMAFD1000	EVOS® FL Auto Imaging System	EA
AMAFD2000	EVOS™ FL Auto 2 Imaging System	EA



Only available for order through legal entity Life Technologies Australia Pty Ltd

Flow Cytometer

Designed with non-clogging fluidics and exceptionally high sample flow rates, the Attune NxT flow cytometer offers rapid, accurate analysis for a broad range of samples—all in a grow-with-you system that can be customized with up to 4 lasers and 14 colors.

Find out more at thermofisher.com/attune

Item	Description	Unit
A24864	Attune NxT Acoustic Focusing Cytometer (Blue laser)	EA
A24863	Attune NxT Acoustic Focusing Cytometer (Blue, red laser)	EA
A24862	Attune NxT Acoustic Focusing Cytometer (Blue, violet laser)	EA
A24861	Attune NxT Acoustic Focusing Cytometer (Blue, yellow laser)	EA
A24860	Attune NxT Acoustic Focusing Cytometer (Blue, red, violet laser)	EA
A24859	Attune NxT Acoustic Focusing Cytometer (Blue, violet, yellow laser)	EA
A24858	Attune NxT Acoustic Focusing Cytometer (Blue, red, yellow, violet laser)	EA



Only available for order through legal entity Life Technologies Australia Pty Ltd

Cellular & Molecular Tools for Wine Analysis

Nucleic Acid Purification

Reap high yields of quality nucleic acids from a variety of starting materials with the Thermo Scientific™ KingFisher™ Duo Prime purification system. This compact instrument is ideal for laboratories requiring nucleic acid purification, and is easily modified for use in proteomics and cell isolation. It provides excellent reproducibility with minimal hands-on time, and has a small footprint in the lab. New features, including a UV lamp for easy decontamination and optional bar-code reader enhance functionality.



Item	Description	Unit
5400110	KingFisher Duo Prime Purification Instrument	EA
AM1836	MagMAX 96 Viral RNA Isolation Kit	EA
AM1939	MagMAX Viral RNA Isolation Kit	EA
AMF4300	EVOS® FL Imaging System	EA
4462359	MagMAX Pathogen RNA/DNA Kit	EA
A25597	MagMAX DNA Multi Sample Ultra Kit	EA
AM1840	MagMax Total Nucleic Acid Isolation Kit	EA

Only available for order through legal entity Life Technologies Australia Pty Ltd

Real-Time PCR Instruments & Assays

Every lab is unique. That's why you deserve a qPCR platform that fits your needs. Perhaps you're looking for simplicity on a budget, or reliable results from limited samples. Maybe your research requires high throughput for maximum productivity, or absolute answers to take your work to the next level. Whatever you need, there's a Quantstudio qPCR system that's just right for you.

Find out more at thermofisher.com/quantstudioqpcrfamily

Item	Description	Unit
A28136	QuantStudio® 3 Real-Time PCR System, 96-well, 0.1mL	EA
A28137	QuantStudio® 3 Real-Time PCR System, 96-well, 0.2mL	EA
A28138	QuantStudio® 5 Real-Time PCR System, 96-well, 0.1mL	EA
A28139	QuantStudio® 5 Real-Time PCR System, 96-well, 0.2mL	EA
A28140	QuantStudio® 5 Real-Time PCR System, 384-well	EA
4396838	TaqMan® Environmental Master Mix 2.0	200 reactions
4444432	TaqMan® Fast Virus 1-Step Master Mix	1 x 1mL
4415440	MeltDoctor™ HRM Master Mix	50mL

Only available for order through legal entity Life Technologies Australia Pty Ltd





Order Placement:

For customer service, call 1300-735-292

To fax an order, use 1800-067-639

To email an order, ordersau@thermofisher.com

Contact Us:

For customer service, call 1300-735-292

For service and calibration, call 1300-736-767

To order online: thermofisher.com.au/godirect

Find out more at thermofisher.com.au

ThermoFisher
S C I E N T I F I C